

Anti-GEMC1 antibody ab112571

★★★★☆ [2 Abreviews](#) [1 图像](#)

概述

产品名称	Anti-GEMC1抗体
描述	兔多克隆抗体to GEMC1
宿主	Rabbit
经测试应用	适用于: WB
种属反应性	与反应: Mouse
免疫原	Synthetic peptide corresponding to Human GEMC1 (internal sequence). A 19 amino acid peptide corresponding to an internal sequence of Human GEMC1 (NP_001140158). Database link: A6NCL1
阳性对照	Mouse Heart tissue lysate.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	pH: 7.2 Preservative: 0.02% Sodium azide Constituent: 99% PBS
纯度	Immunogen affinity purified
纯化说明	ab112571 is purified by affinity chromatography purified via peptide column.
克隆	多克隆
同种型	IgG

应用

The Abpromise guarantee

Abpromise™承诺保证使用ab112571于以下的经测试应用

“应用说明”部分下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。


应用	Ab评论	说明
WB		Use a concentration of 1 - 2 µg/ml. Predicted molecular weight: 38 kDa.

靶标

相关性 Regulator of DNA replication. Promotes initiation of chromosomal DNA replication by mediating TOPBP1- and CDK2-dependent recruitment of CDC45L onto replication origins.

细胞定位 Nuclear

图片



All lanes : Anti-GEMC1 antibody (ab112571) at 1 µg/ml

Lane 1 : Mouse heart tissue lysate

Lane 2 : Mouse heart tissue lysate with blocking peptide

Lysates/proteins at 15 µg per lane.

Predicted band size: 38 kDa

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors