

# Anti-GM-CSF antibody ab9667

**2 References**   **1 图像**

### 概述

产品名称	Anti-GM-CSF抗体
描述	兔多克隆抗体to GM-CSF
宿主	Rabbit
经测试应用	适用于: Sandwich ELISA
种属反应性	与反应: Recombinant fragment
免疫原	Recombinant full length protein corresponding to Human GM-CSF . Highly pure (>98%) recombinant hGM-CSF (human Granulocyte Macrophage Colony Stimulating Factor) Database link: <b>P04141</b>
阳性对照	Recombinant human GM-CSF protein ( <b>ab9668</b> ) can be used as a positive control in WB.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

### 性能

形式	Lyophilized:Reconstitute with 200µl of sterile water. Please note that if you receive this product in liquid form it has already been reconstituted as described and no further reconstitution is necessary.
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C.
存储溶液	No preservative, sterile filtered
无载体	是
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

### 应用

The Abpromise guarantee **Abpromise™**承诺保证使用ab9667于以下的经测试应用


“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
Sandwich ELISA		Use a concentration of 0.5 - 2 µg/ml. Allows the detection of 0.2 - 0.4 ng/well of recombinant hGM-CSF (using 100µl/well antibody solution). Can be paired for Sandwich ELISA with <a href="#">ab84259</a> . Can be used as the capture antibody when paired with recommended pair.

靶标

功能	Cytokine that stimulates the growth and differentiation of hematopoietic precursor cells from various lineages, including granulocytes, macrophages, eosinophils and erythrocytes.
序列相似性	Belongs to the GM-CSF family.
细胞定位	Secreted.

图片



Sandwich ELISA of ab9667. Allows the detection of 0.2 - 0.4 ng/well of recombinant hGM-CSF (using 100µl/well antibody solution). Can be paired for Sandwich ELISA with [ab84259](#). Can be used as the capture antibody when paired with recommended pair.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish

- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

#### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors