

### Anti-GPR91 antibody ab75105

★★★★☆ [1 Abreviews](#) [1 References](#) [1 图像](#)

#### 概述

产品名称	Anti-GPR91抗体
描述	兔多克隆抗体to GPR91
宿主	Rabbit
经测试应用	适用于: WB
种属反应性	与反应: Human 预测可用于: Mouse, Rat 
免疫原	Synthetic peptide (Human)(Peptide available as <a href="#">ab106670</a> .)
阳性对照	Extracts from MCF-7 and HUVEC cells.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

#### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at -20°C. Stable for 12 months at -20°C.
存储溶液	pH: 7.40 Preservative: 0.02% Sodium azide Constituents: 50% Glycerol (glycerin, glycerine), 0.87% Sodium chloride, PBS
纯度	Without Mg <sup>2+</sup> and Ca <sup>2+</sup>
克隆	Immunogen affinity purified
同种型	多克隆 IgG

#### 应用

## The Abpromise guarantee

**Abpromise™**承诺保证使用ab75105于以下的经测试应用


“应用说明”部分下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB	★★★★☆ (1)	1/500 - 1/1000. Detects a band of approximately 38 kDa (predicted molecular weight: 38 kDa).

## 靶标

功能	Receptor for succinate.
组织特异性	Expressed specifically in kidney.
序列相似性	Belongs to the G-protein coupled receptor 1 family.
细胞定位	Cell membrane.

## 图片



**All lanes** : Anti-GPR91 antibody (ab75105) at 1/500 dilution

**Lane 1** : Extracts from HUVEC cells

**Lane 2** : Extracts from MCF-7 cells

**Lane 3** : Extracts from MCF-7 cells with immunising peptide at 5  $\mu$ g

Lysates/proteins at 5  $\mu$ g per lane.

**Predicted band size:** 38 kDa

**Observed band size:** 38 kDa

**Additional bands at:** 100 kDa, 30 kDa. We are unsure as to the identity of these extra bands.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you

- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

#### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors