

Anti-Granzyme A antibody [GA6] ab10870

★★★★★ [1 Abreviews](#) [4 References](#)

概述

产品名称	Anti-Granzyme A抗体[GA6]
描述	小鼠单克隆抗体[GA6] to Granzyme A
宿主	Mouse
经测试应用	适用于: IHC-P, IHC-FoFr, WB
种属反应性	与反应: Human
免疫原	Recombinant full length protein corresponding to Human Granzyme A. Database link: P12544
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: 0.09% Sodium azide Constituent: PBS
纯度	Protein A purified
纯化说明	Purified IgG prepared by affinity chromatography from tissue culture supernatant.
克隆	单克隆
克隆编号	GA6
骨髓瘤	Sp2/0
同种型	IgG1

The Abpromise guarantee **Abpromise™**承诺保证使用ab10870于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
IHC-P	★★★★★ (1)	Use at an assay dependent concentration.
IHC-FoFr		Use at an assay dependent concentration. PubMed: 19827750
WB		Use at an assay dependent concentration. Predicted molecular weight: 29 kDa. Recognises a 60kD disulphide-linked homodimeric protein.

靶标

功能	This enzyme is necessary for target cell lysis in cell-mediated immune responses. It cleaves after Lys or Arg. Cleaves APEX1 after 'Lys-31' and destroys its oxidative repair activity. Involved in apoptosis.
序列相似性	Belongs to the peptidase S1 family. Granzyme subfamily. Contains 1 peptidase S1 domain.
细胞定位	Secreted. Cytoplasmic granule.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors