

### Anti-HAI-1 antibody [MM0322-3X28] ab89404

#### 1 图像

#### 概述

产品名称	Anti-HAI-1抗体[MM0322-3X28]
描述	小鼠单克隆抗体[MM0322-3X28] to HAI-1
宿主	Mouse
经测试应用	适用于: WB
种属反应性	与反应: Human
免疫原	Recombinant fragment corresponding to Human HAI-1.
常规说明	<p>Reconstituted antibody should be aliquoted and can be stored frozen at &lt; -20°C for at least six months without detectable loss of activity. Avoid repeated freeze-thaw cycles.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

#### 性能

形式	Lyophilized:Reconstitute the antibody with 200 µl sterile PBS to give a final concentration of 0.5mg/ml.
存放说明	Shipped at 4°C. Store at -80°C.
存储溶液	Constituent: PBS
纯度	Protein G purified
纯化说明	The IgG fraction of culture supernatant was purified by Protein G affinity chromatography and filtered through a 0.2 µm filter.
克隆	单克隆
克隆编号	MM0322-3X28
同种型	IgG1

## The Abpromise guarantee Abpromise™承诺保证使用ab89404于以下的经测试应用


“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB		1/500 - 1/1000. Predicted molecular weight: 58 kDa.

## 靶标

功能	Inhibitor of HGF activator. Also acts as an inhibitor of matriptase (ST14).
序列相似性	Contains 2 BPTI/Kunitz inhibitor domains. Contains 1 LDL-receptor class A domain. Contains 1 MANSC domain.
结构域	This inhibitor contains two inhibitory domains.
细胞定位	Secreted.

## 图片



Anti-HAI-1 antibody [MM0322-3X28] (ab89404) at 1/1000 dilution  
(Mouse anti human HAI-1) + Lysates of human MCF-7 cells

**Predicted band size:** 58 kDa

**Observed band size:** 65 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish

- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

#### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors