

Anti-Hairy antibody [1/24] ab20165

★ ★ ★ ★ ★ [1 Abreviews](#) [3 References](#)

概述

产品名称	Anti-Hairy抗体[1/24]
描述	小鼠单克隆抗体[1/24] to Hairy
宿主	Mouse
经测试应用	适用于: ChIP, IHC, WB
种属反应性	与反应: Drosophila melanogaster
免疫原	Other Immunogen Type corresponding to Hairy.
常规说明	<p>This antibody stains Drosophila embryos.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: 0.02% Sodium azide Constituent: 99.98% PBS
纯度	Protein A/G purified
Primary antibody说明	This antibody stains Drosophila embryos
克隆	单克隆
克隆编号	1/24
骨髓瘤	Sp2/0-Ag14
同种型	IgG1
轻链类型	unknown

应用

The Abpromise guarantee **Abpromise™**承诺保证使用ab20165于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ChIP		Use at an assay dependent concentration.
IHC		Use at an assay dependent concentration.
WB		Use at an assay dependent concentration. Predicted molecular weight: 37 kDa.

靶标

相关性	Pair rule protein that regulates embryonic segmentation and adult bristle patterning. Transcriptional repressor of genes that require a bHLH protein for their transcription (e.g. the Fushi tarazu gene. Transcription repression requires formation of a complex with a corepressor
细胞定位	Nuclear

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors