

Anti-Hemopexin antibody [04] ab27710

概述

产品名称	Anti-Hemopexin抗体[04]
描述	小鼠单克隆抗体[04] to Hemopexin
宿主	Mouse
经测试应用	适用于: ELISA, WB
种属反应性	与反应: Human
免疫原	Full length native protein: Hemopexin purified from human plasma. The immunogen has been adsorbed onto Al(OH) ₃ prior to injection in mice.
常规说明	<p>The fusion partner for this antibody is SP2mIL6.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at +4°C. Avoid freeze / thaw cycle.
存储溶液	<p>pH: 7.40</p> <p>Preservative: 0.097% Sodium azide</p> <p>Constituents: 0.0268% PBS, 2.9% Sodium chloride</p>
纯度	Protein A purified
Primary antibody说明	The fusion partner for this antibody is SP2mIL6.
克隆	单克隆
克隆编号	04
同种型	IgG2b
轻链类型	kappa

应用

The Abpromise guarantee **Abpromise™**承诺保证使用ab27710于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ELISA		Use at an assay dependent concentration.
WB		Use at an assay dependent concentration.

靶标

功能	Binds heme and transports it to the liver for breakdown and iron recovery, after which the free hemopexin returns to the circulation.
组织特异性	Expressed by the liver and secreted in plasma.
序列相似性	Belongs to the hemopexin family. Contains 5 hemopexin-like domains.
翻译后修饰	N- and O-glycosylated. O-glycosylated with core 1 or possibly core 8 glycans.
细胞定位	Secreted.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors