

### Anti-Hemopexin antibody [32] ab27711

#### 1 图像

#### 概述

产品名称	Anti-Hemopexin抗体[32]
描述	小鼠单克隆抗体[32] to Hemopexin
宿主	Mouse
经测试应用	适用于: Sandwich ELISA
种属反应性	与反应: Human
免疫原	Full length native protein: Hemopexin purified from human plasma. Immunogen was adsorbed onto Al(OH) <sub>3</sub> prior to injection in mice.
常规说明	<p>The fusion partner for this antibody is SP2mIL6.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

#### 性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	<p>pH: 7.40</p> <p>Preservative: 0.097% Sodium azide</p> <p>Constituents: 0.0268% PBS, 2.9% Sodium chloride</p>
纯度	Protein A purified
Primary antibody说明	The fusion partner for this antibody is SP2mIL6.
克隆	单克隆
克隆编号	32
同种型	IgG1
轻链类型	kappa

## 应用

**The Abpromise guarantee**      **Abpromise™**承诺保证使用ab27711于以下的经测试应用


“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
<b>Sandwich ELISA</b>		Use a concentration of 5 µg/ml. Can be paired for Sandwich ELISA with <b><u>Rabbit polyclonal to Hemopexin (ab48135)</u></b> . For sandwich ELISA, use this antibody as Capture at 5 µg/ml with <b><u>Rabbit polyclonal to Hemopexin (ab48135)</u></b> as Detection.

## 靶标

<b>功能</b>	Binds heme and transports it to the liver for breakdown and iron recovery, after which the free hemopexin returns to the circulation.
<b>组织特异性</b>	Expressed by the liver and secreted in plasma.
<b>序列相似性</b>	Belongs to the hemopexin family. Contains 5 hemopexin-like domains.
<b>翻译后修饰</b>	N- and O-glycosylated. O-glycosylated with core 1 or possibly core 8 glycans.
<b>细胞定位</b>	Secreted.

## 图片



Standard Curve for Hemopexin (Analyte: **Hemopexin protein (ab98988)**); dilution range 1pg/ml to 1µg/ml using Capture Antibody **Mouse monoclonal [32] to Hemopexin (ab27711)** at 5µg/ml and Detector Antibody **Rabbit polyclonal to Hemopexin (ab48135)** at 0.5µg/ml.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet

- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors