

Anti-Hepatitis A Virus antibody ab68579

概述

产品名称	Anti-Hepatitis A Virus抗体
描述	山羊多克隆抗体to Hepatitis A Virus
宿主	Goat
特异性	This antibody reacts with Hepatitis A virions in viral and infected tissue samples.
经测试应用	适用于: ELISA
种属反应性	与反应: Hepatitis A virus
免疫原	Tissue, cells or virus corresponding to Hepatitis A Virus. Hepatitis A virus strain HM175
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	<p>pH: 7.2</p> <p>Preservative: 0.1% Sodium azide</p> <p>Constituent: 0.0268% PBS</p>
纯度	Affinity purified
克隆	多克隆
同种型	IgG

应用

The Abpromise guarantee **Abpromise™** 承诺保证使用ab68579于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ELISA		Use at an assay dependent concentration.

靶标

相关性	Hepatitis A Virus (HAV) is a 27nm nonenveloped, spherical, positive stranded RNA virus, classified within the genus hepatovirus of the picornavirus family and is among the smallest and structurally simplest of the RNA animal viruses. A single large polyprotein is expressed from a large open reading frame extending through most of the genomic RNA. This polyprotein is subsequently cleaved by a viral protease (3Cpro) to form three (possibly four) capsid proteins and several nonstructural proteins. HAV genomic replication occurs exclusively in the cytoplasm of the infected hepatocyte by a mechanism involving an RNA-dependent RNA polymerase.
-----	--

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors