

Anti-Hexokinase II antibody [3D3] ab104836

★★★★★ [8 Abreviews](#) [66 References](#) [6 图像](#)

概述

产品名称	Anti-Hexokinase II抗体[3D3]
描述	小鼠单克隆抗体[3D3] to Hexokinase II
宿主	Mouse
经测试应用	适用于: WB, IHC-P, Flow Cyt
种属反应性	与反应: Human
免疫原	Recombinant fragment corresponding to Human Hexokinase II.
阳性对照	WB: Jurkat, HeLa and HEK293 cell lysates. IHC-P: Human esophagus and lung cancer tissues. Flow Cyt: K562 cells.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: 0.03% Sodium azide Constituent: PBS
纯度	Protein G purified
纯化说明	Purified from tissue culture supernatant.
克隆	单克隆
克隆编号	3D3
同种型	IgG1

应用

The Abpromise guarantee


Abpromise™承诺保证使用ab104836于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB	★★★★★ (4)	1/500 - 1/2000. Predicted molecular weight: 102 kDa.
IHC-P	★★★★☆ (3)	1/200 - 1/1000.
Flow Cyt		Use at an assay dependent concentration. ab170190 - Mouse monoclonal IgG1, is suitable for use as an isotype control with this antibody.

靶标	
组织特异性	Predominant hexokinase isozyme expressed in insulin-responsive tissues such as skeletal muscle.
通路	Carbohydrate metabolism; hexose metabolism.
序列相似性	Belongs to the hexokinase family. Contains 2 hexokinase domains.
结构域	The N- and C-terminal halves of this hexokinase show extensive sequence similarity to each other. The catalytic activity is associated with the C-terminus while regulatory function is associated with the N-terminus. Each domain can bind a single glucose and Gluc-6-P molecule.
细胞定位	Mitochondrion outer membrane. Its hydrophobic N-terminal sequence may be involved in membrane binding.

图片



Western blot - Anti-Hexokinase II antibody [3D3] (ab104836)

All lanes : Anti-Hexokinase II antibody [3D3] (ab104836) at 1/500 dilution

Lane 1 : Wild-type HEK-293T cell lysate

Lane 2 : HK2 CRISPR/Cas9 edited HEK-293T cell lysate

Lane 3 : HeLa cell lysate

Lane 4 : A549 cell lysate

Lysates/proteins at 20 µg per lane.


Performed under reducing conditions.

Predicted band size: 102 kDa

Observed band size: 102 kDa

Lanes 1 - 4: Merged signal (red and green). Green - ab104836 observed at 102 kDa. Red - loading control **ab181602** (Rabbit Anti-GAPDH antibody [EPR16891]) observed at 37kDa.

ab104836 was shown to react with Hexokinase II in wild-type HEK-293T cells in western blot. The band observed in HK2 CRISPR/Cas9 edited sample lane below 102 kDa may represent truncated forms and cleaved fragments. This has not been investigated further. Wild-type and HK2 CRISPR/Cas9 edited HEK-293T cell lysates were subjected to SDS-PAGE. Membranes were blocked in 3% milk in TBS-T (0.1% Tween[®]) before incubation with ab104836 and **ab181602** (Rabbit Anti-GAPDH antibody [EPR16891]) overnight at 4°C at a 1 in 500 Dilution and a 1 in 20000 dilution respectively. Blots were incubated with Goat anti-Mouse IgG H&L (IRDye[®] 800CW) preabsorbed (**ab216772**) and Goat anti-Rabbit IgG H&L (IRDye[®] 680RD) preabsorbed (**ab216777**) secondary antibodies at 1 in 20000 dilution for 1 hour at room temperature before imaging.




Immunohistochemical analysis of human skin tissue, staining Hexokinase II with ab104836.

Tissue was fixed with paraformaldehyde and blocked with 10% BSA for 2 hours at room temperature; antigen retrieval was by heat mediation in citrate buffer. Samples were incubated with primary antibody (1/250 in 5% NGS + 10% BSA in PBS) for 16 hours at 4°C. An HRP-conjugated goat anti-mouse polyclonal IgG (1/250) was used as the secondary antibody.


Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Hexokinase II antibody [3D3] (ab104836)

This image is courtesy of an Abreview submitted by Ms Hadas Pahima



Flow Cytometry - Anti-Hexokinase II antibody [3D3]
(ab104836)

ab104836 staining Hexokinase II in K562 cells by Immunofluorescence (green) and negative control (purple).



Western blot - Anti-Hexokinase II antibody [3D3]
(ab104836)


All lanes : Anti-Hexokinase II antibody [3D3] (ab104836) at 1/500 dilution

Lane 1 : Jurkat cell lysate

Lane 2 : HeLa cell lysate


Lane 3 : HEK293 cell lysate

Predicted band size: 102 kDa



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Hexokinase II antibody [3D3] (ab104836)

ab104836 at 1/200 dilution staining Hexokinase II in human esophagus cancer by Immunohistochemistry, Paraffin-embedded tissue with DAB staining.



ab104836 at 1/200 dilution staining Hexokinase II in human lung cancer by Immunohistochemistry, Paraffin-embedded tissue with DAB staining.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Hexokinase II antibody [3D3] (ab104836)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors