

# Anti-Hirudin antibody [2D7] ab17198

★★★★★ [1 Abreviews](#) [1 References](#)

## 概述

产品名称	Anti-Hirudin抗体[2D7]
描述	小鼠单克隆抗体[2D7] to Hirudin
宿主	Mouse
特异性	The affinity (Kd) of this antibody for rHV2K47 is 4.5x10E-7 M
经测试应用	适用于: ELISA, WB
种属反应性	与反应: Leech
免疫原	Recombinant full length protein corresponding to Leech Hirudin.
表位	The epitope of this antibody is located at the C-terminal end of hirudin.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

## 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	<p>pH: 7.40</p> <p>Preservative: 0.097% Sodium azide</p> <p>Constituents: 0.0268% PBS, 2.9% Sodium chloride</p>
纯度	Protein G purified
克隆	单克隆
克隆编号	2D7
骨髓瘤	x63-Ag8.653
同种型	IgG1
轻链类型	kappa

## 应用

**The Abpromise guarantee**      **Abpromise™**承诺保证使用ab17198于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ELISA		1/2000.
WB	★★★★★ (1)	Use at an assay dependent concentration. Detects a band of approximately 7 kDa. ab17198 reacts with hirudin in both reduced and unreduced forms.

## 靶标

相关性	Hirudin is a highly effective and very specific anticoagulant which specifically inhibits formation of fibrin and activation of factors V, VIII and XIII by forming irreversible, noncovalent 1:1 complexes with a-thrombin. Hirudin is a 65 amino acid polypeptide, 7 kDa protein naturally originating from the medical leech Hirudo medicinalis. It is the most potent thrombin inhibitor known to date.
细胞定位	Secreted

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors