

# Anti-HOXA7 antibody [743C1a] ab51235

**2 References**   **1 图像**

### 概述

产品名称	Anti-HOXA7抗体[743C1a]
描述	小鼠单克隆抗体[743C1a] to HOXA7
宿主	Mouse
经测试应用	适用于: WB
种属反应性	与反应: Recombinant fragment
免疫原	Recombinant fragment corresponding to Human HOXA7.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

### 性能

形式	Liquid
存放说明	Store at +4°C short term (1-2 weeks) (add antibacterial agent). Aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	<p>pH: 7.40</p> <p>Preservative: 0.05% Sodium azide</p> <p>Constituents: 1% BSA, PBS</p>
纯度	Protein G purified
纯化说明	ab51235 was purified from culture supernatant of the hybridoma, cultured in a medium containing bovine IgG-depleted (approximately 95%) fetal bovine serum.
克隆	单克隆
克隆编号	743C1a
同种型	IgG1

### 应用

## The Abpromise guarantee

**Abpromise™** 承诺保证使用ab51235于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB		Use at an assay dependent concentration. Predicted molecular weight: 25 kDa.

## 靶标

功能	Sequence-specific transcription factor which is part of a developmental regulatory system that provides cells with specific positional identities on the anterior-posterior axis.
序列相似性	Belongs to the Antp homeobox family. Contains 1 homeobox DNA-binding domain.
细胞定位	Nucleus.

## 图片



Anti-HOXA7 antibody [743C1a] (ab51235)

**Predicted band size:** 25 kDa

**Observed band size:** 31 kDa

Western blot analysis of immunised recombinant protein

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors