

# Anti-Hrk antibody ab45419

★★★★☆ [1 Abreviews](#) [5 References](#) [1 图像](#)

### 概述

|       |   |
|-------|---|
| 产品名称  | Anti-Hrk抗体  |
| 描述    | 兔多克隆抗体to Hrk  |
| 宿主    | Rabbit  |
| 经测试应用 | 适用于: ICC/IF, WB   |
| 种属反应性 | 与反应: Mouse, Human   |
| 免疫原   | Synthetic peptide of 15 amino acids in length from near the center of human Hrk   |
| 阳性对照  | Mouse pancreas tissue lysate.   |
| 常规说明  | <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p> |

### 性能

|      |   |
|------|---|
| 形式   | Liquid  |
| 存放说明 | Shipped at 4°C. Store at +4°C.                                  |
| 存储溶液 | pH: 7.2<br>Preservative: 0.02% Sodium azide<br>Constituent: PBS |
| 纯度   | Immunogen affinity purified                                     |
| 克隆   | 多克隆   |
| 同种型  | IgG   |

### 应用

**The Abpromise guarantee**      **Abpromise™**承诺保证使用ab45419于以下的经测试应用


“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

| 应用     | Ab评论      | 说明   |
|--------|-----------|--|
| ICC/IF |           | Use at an assay dependent concentration.   |
| WB     | ★★★★☆ (1) | Use a concentration of 2.5 - 5 µg/ml. Detects a band of approximately 13 kDa (predicted molecular weight: 10 kDa). |

## 靶标

|      |  |
|------|--|
| 功能   | Promotes apoptosis.  |
| 结构域  | The BH3 motif is required for the induction of cell death. |
| 细胞定位 | Membrane. Mitochondrion.                                   |

## 图片



**Lane 1** : Anti-Hrk antibody (ab45419) at 2.5 µg/ml

**Lane 2** : Anti-Hrk antibody (ab45419) at 5 µg/ml

**All lanes** : Mouse pancreas tissue lysate

**Predicted band size:** 10 kDa

**Observed band size:** 13 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors