

# HRP Anti-Bovine Serum Albumin antibody ab7637

★★★★★ [3 Abreviews](#) [1 References](#)

### 概述

---

产品名称	HRP Anti-Bovine Serum Albumin抗体
描述	HRP兔多克隆抗体to Bovine Serum Albumin
宿主	Rabbit
偶联物	HRP
经测试应用	适用于: ELISA, IHC-Fr, Dot blot, Immunomicroscopy, WB
种属反应性	与反应: Cow
免疫原	Full length native Albumin (purified).
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

### 性能

---

形式	Liquid
存放说明	Store at +4°C. Store In the Dark.
存储溶液	Preservative: 0.01% Gentamicin sulphate Constituents: 1% Polyethylene glycol, 0.42% Potassium phosphate, 0.87% Sodium chloride
纯度	IgG fraction
纯化说明	This product is an IgG fraction antibody purified from monospecific antiserum by a multi-step process which includes delipidation, salt fractionation and ion exchange chromatography followed by extensive dialysis.
克隆	多克隆
同种型	IgG

### 应用

---

## The Abpromise guarantee **Abpromise™**承诺保证使用ab7637于以下的经测试应用

“应用说明”部分下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ELISA		Use at an assay dependent concentration.
IHC-Fr		Use at an assay dependent concentration.
Dot blot		Use at an assay dependent concentration.
Immunomicroscopy		Use at an assay dependent concentration.
WB	★★★★★ (3)	Use at an assay dependent concentration.

## 靶标

---

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors