

### HRP Anti-Cathelicidin/CLP antibody [H7] ab187933

#### 概述

---

<b>产品名称</b>	HRP Anti-Cathelicidin/CLP抗体[H7]
<b>描述</b>	HRP小鼠单克隆抗体[H7] to Cathelicidin/CLP
<b>宿主</b>	Mouse
<b>偶联物</b>	HRP
<b>经测试应用</b>	<b>适用于:</b> WB, Flow Cyt, ICC/IF, IP
<b>种属反应性</b>	<b>与反应:</b> Human
<b>免疫原</b>	Full length native protein (purified) corresponding to Human Cathelicidin/CLP. (Partially purified human neutrophil membrane proteins (heparin Ultrogellectin affinity purification fractions in octyl glucoside)).
<b>表位</b>	ab187933 recognises amino acid sequences 115-NQARGSF-121 in the pro-domain of cathelicidin.
<b>常规说明</b>	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

#### 性能

---

<b>形式</b>	Liquid
<b>存放说明</b>	Shipped at 4°C. Store at +4°C.
<b>存储溶液</b>	Preservative: 0.01% Thimerosal (merthiolate) Constituent: 99% PBS
<b>纯度</b>	Protein G purified
<b>克隆</b>	单克隆
<b>克隆编号</b>	H7
<b>骨髓瘤</b>	P3U1
<b>同种型</b>	IgG1

## 应用

### The Abpromise guarantee **Abpromise™**承诺保证使用ab187933于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB		1/1000. Predicted molecular weight: 19 kDa.
Flow Cyt		Use at an assay dependent concentration.
ICC/IF		Use at an assay dependent concentration.
IP		Use at an assay dependent concentration.

## 靶标

功能	Binds to bacterial lipopolysaccharides (LPS), has antibacterial activity.
组织特异性	Expressed in bone marrow and testis and neutrophils.
序列相似性	Belongs to the cathelicidin family.
翻译后修饰	The N-terminus is blocked.
细胞定位	Secreted.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors