

HRP Anti-Diphtheria Toxin antibody ab34420

概述

产品名称	HRP Anti-Diphtheria Toxin抗体
描述	HRP山羊多克隆抗体to Diphtheria Toxin
宿主	Goat
偶联物	HRP
特异性	Reacts with purified diphtheria toxin/toxoid.
经测试应用	适用于: ELISA, WB
种属反应性	与反应: Corynebacterium diphtheriae
免疫原	Full length protein corresponding to Diphtheria Toxin. Full length purified diphtheria toxoid protein from Corynebacterium diphtheriae. Database link: Q6NK15
常规说明	<p>Highly purified preparation of horseradish peroxidase. Free enzyme is absent.</p> <p>Storage in frost-free freezers is not recommended.</p> <p>This product should be stored undiluted. Should this product contain a precipitate we recommend microcentrifugation before use.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C. Avoid freeze / thaw cycle. Please see notes section.
存储溶液	<p>pH: 7.20</p> <p>Preservative: 0.002% Thimerosal (merthiolate)</p> <p>Constituents: 1% BSA, PBS</p>
纯化说明	Purified IgG
克隆	多克隆

应用

The Abpromise guarantee

Abpromise™承诺保证使用ab34420于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ELISA		1/200 - 1/1000.
WB		1/20 - 1/200. Predicted molecular weight: 58 kDa.

靶标

相关性

Corynebacterium diphtheriae is a gram-positive, nonmotile bacteria found in soil and animal feces. C. diphtheriae infect the epithelial cells of the upper respiratory tract from where they produce and secrete a potent toxin. This toxin is absorbed and disseminated through lymph channels and blood to the susceptible tissues of the body.

细胞定位

Secreted

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors