

HRP Anti-E tag antibody ab19400

[5 References](#) [1 图像](#)

概述

产品名称	HRP Anti-E tag抗体
描述	HRP山羊多克隆抗体to E tag
宿主	Goat
偶联物	HRP
经测试应用	适用于: ICC, ELISA, WB
种属反应性	与反应: Species independent
免疫原	Synthetic peptide: GAPVPYPDPLEPR (E tag) conjugated to KLH.

 [Run BLAST with](#)

 [Run BLAST with](#)

常规说明

The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.

If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	pH: 6.8 Preservative: 0.01% Thimerosal (merthiolate) Constituents: 1.19% HEPES, 0.58% Sodium chloride, 0.2% BSA
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

The Abpromise guarantee

Abpromise™ 承诺保证使用 ab19400 于以下的经测试应用

“应用说明”部分下显示的仅为推荐的起始稀释度；实际最佳的稀释度/浓度应由使用者检定。


应用	Ab评论	说明
ICC		1/200 - 1/500.
ELISA		1/10000 - 1/100000.
WB		1/1000 - 1/10000. Predicted molecular weight: 1 kDa. Colorimetric: 1/1000 - 1/10000. Chemiluminescence: 1/1000 - 1/30000.

靶标

相关性

E tag is a peptide epitope tag, sequence GAPVPYPDPLEPR.

图片



All lanes : HRP Anti-E tag antibody (ab19400) at 0.04 µg/ml

Lane 1 : E. coli whole cell lysate expressing a multi-tag fusion protein at 0.2 µg

Lane 2 : E. coli whole cell lysate expressing a multi-tag fusion protein at 0.1 µg

Lane 3 : E. coli whole cell lysate expressing a multi-tag fusion protein at 0.05 µg

Predicted band size: 1 kDa

Exposure time: 3 seconds

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors