

HRP Anti-EGF antibody [KT2] ab106489

1 图像

概述	
产品名称	HRP Anti-EGF 抗体[KT2]
描述	HRP 小鼠单克隆抗体[KT2] to EGF
宿主	Mouse
偶联物	HRP
经测试应用	适用于: Sandwich ELISA
种属反应性	与反应: Human
免疫原	Recombinant mature full length protein (Human) (53aa)
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>
性能	
形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	Preservative: 0.01% Thimerosal (merthiolate) Constituent: PBS
纯度	Protein A purified
克隆	单克隆
克隆编号	KT2
同种型	IgG2a
应用	

The Abpromise guarantee **Abpromise™** 承诺保证使用ab106489于以下的经测试应用


“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
Sandwich ELISA		1/1000.

靶标

功能	EGF stimulates the growth of various epidermal and epithelial tissues in vivo and in vitro and of some fibroblasts in cell culture. Magnesiotropic hormone that stimulates magnesium reabsorption in the renal distal convoluted tubule via engagement of EGFR and activation of the magnesium channel TRPM6. Can induce neurite outgrowth in motoneurons of the pond snail <i>Lymnaea stagnalis</i> in vitro (PubMed:10964941).
组织特异性	Expressed in kidney, salivary gland, cerebrum and prostate.
疾病相关	Hypomagnesemia 4
序列相似性	Contains 9 EGF-like domains. Contains 9 LDL-receptor class B repeats.
翻译后修饰	O-glycosylated with core 1-like and core 2-like glycans. It is uncertain if Ser-954 or Thr-955 is O-glycosylated. The modification here shows glycan heterogeneity: HexHexNAc (major) and Hex2HexNAc2 (minor).
细胞定位	Membrane.

图片



Sandwich ELISA - HRP Anti-EGF antibody [KT2]
(ab106489)

Sandwich ELISA using a mouse monoclonal antibody to EGF coated plate and ab106489 at 1/1000 dilution.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish

- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors