

HRP Anti-HIV1 gp41 antibody ab68622

1 References

概述

产品名称	HRP Anti-HIV1 gp41抗体
描述	HRP山羊多克隆抗体to HIV1 gp41
宿主	Goat
偶联物	HRP
特异性	This antibody reacts with HIV1 gp41 in viral and infected tissue samples.
经测试应用	适用于: WB, ELISA
种属反应性	与反应: Human immunodeficiency virus
免疫原	Recombinant fragment, corresponding to the ectodomain of gp41-glycosylated (a.a. residues 546-682) from strain HxB2.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	Preservative: 0.002% Thimerosal (merthiolate) Constituents: 1% BSA, PBS
纯度	Affinity purified
克隆	多克隆
同种型	IgG

应用

Abpromise™ 承诺保证使用ab68622于以下的经测试应用

The Abpromise guarantee

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB		Use at an assay dependent concentration. Predicted molecular weight: 41 kDa.
ELISA		1/200 - 1/1000.

靶标

相关性

gp41/120 is the major HIV protein associated with the HIV envelope. It functions as the viral antireceptor or attachment protein. gp41 (or TM) traverses the envelope, whereas gp120 is present on the outer surface and is noncovalently attached to gp41. The precursor of gp120/41 (gp160) is synthesized in the endoplasmic reticulum and is transported via the golgi body to the cell surface. Upon activation of the envelope glycoprotein (gp120/41) by cellular receptors, gp41 undergoes conformational changes that mediate fusion of the viral and cellular membranes.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors