

HRP Anti-IL-11 antibody [KT62] ab106492

1 图像

概述

产品名称	HRP Anti-IL-11抗体[KT62]
描述	HRP小鼠单克隆抗体[KT62] to IL-11
宿主	Mouse
偶联物	HRP
特异性	ab106492 is specific for Human Interleukin 11
经测试应用	适用于: Sandwich ELISA
种属反应性	与反应: Human
免疫原	Recombinant full length protein corresponding to Human IL-11. Sequence: MNCVCRLVLV VLSLWPD TAV APGPPPGPPR VSPDPRAELD STVLLTRSL L ADTRQLAAQL RDKFPADGDH NLDSLPTL AM SAGALGALQL PGVLTRLRAD LLSYLRHVQW LRRAGGSSLK TLEPELGTLQ ARLDRLLRRL QLLMSRLALP QPPPDPPAPP LAPPSSAWGG IRAAHAI

 [Run BLAST with](#)

 [Run BLAST with](#)

常规说明

The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.

If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	Preservative: 0.01% Thimerosal (merthiolate) Constituent: PBS
纯度	Protein A purified

克隆	单克隆
克隆编号	KT62
同种型	IgG1

应用

The Abpromise guarantee **Abpromise™**承诺保证使用ab106492于以下的经测试应用


“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
Sandwich ELISA		1/1000.

靶标

功能	Directly stimulates the proliferation of hematopoietic stem cells and megakaryocyte progenitor cells and induces megakaryocyte maturation resulting in increased platelet production.
序列相似性	Belongs to the IL-6 superfamily.
细胞定位	Secreted.

图片



ab106492 at 1/1000 dilution detecting IL-11 by sandwich ELISA.
The capture antibody is an mouse monoclonal (IgG1) raised against Human IL-11.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish

- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors