

### HRP Anti-IL-17A antibody ab106138

#### 1 图像

#### 概述

产品名称	HRP Anti-IL-17A抗体
描述	HRP兔多克隆抗体to IL-17A
宿主	Rabbit
偶联物	HRP
特异性	In ELISA and other immunoreactive assays, ab106138 will recognize both native and recombinant IL-17 in cell supernatants and certain body fluids.
经测试应用	适用于: WB
种属反应性	与反应: Human
免疫原	Recombinant full length protein corresponding to Human IL-17A.
阳性对照	Recombinant Human IL-17a protein.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

#### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	<p>Preservative: 0.01% Gentamicin sulphate</p> <p>Constituents: 0.42% Potassium phosphate, 0.88% Sodium chloride, 1% BSA</p> <p>BSA is immunoglobulin and protease free. Do NOT add Sodium Azide!</p>
纯度	Protein A purified
纯化说明	This purified antibody was prepared from whole rabbit serum, and has been heated to 56°C for 30 minutes.
克隆	多克隆

同种型

IgG

应用

The Abpromise guarantee

Abpromise™承诺保证使用ab106138于以下的经测试应用


“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB		1/1000 - 1/5000. Predicted molecular weight: 17 kDa.

靶标

功能	Induces stromal cells to produce proinflammatory and hematopoietic cytokines. Enhances the surface expression of the intracellular adhesion molecule-1 (ICAM-1) in fibroblasts.
组织特异性	Restricted to activated memory T-cells.
序列相似性	Belongs to the IL-17 family.
翻译后修饰	Found both in glycosylated and nonglycosylated forms.
细胞定位	Secreted.

图片



HRP Anti-IL-17A antibody (ab106138) at 1/1000 dilution + recombinant Human IL-17A

**Secondary**  
DyLight 649 conjugated anti-Rabbit IgG at 1/20000 dilution

**Predicted band size:** 17 kDa  
**Observed band size:** 17 kDa

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors