

# HRP Anti-Staphylococcus aureus antibody ab68952

### 概述

---

<b>产品名称</b>	HRP抗金黄色葡萄球菌抗体
<b>描述</b>	HRP兔多克隆抗体to金黄色葡萄球菌
<b>宿主</b>	Rabbit
<b>偶联物</b>	HRP
<b>经测试应用</b>	<b>适用于:</b> WB, ELISA
<b>种属反应性</b>	<b>与反应:</b> Staphylococcus aureus
<b>免疫原</b>	Tissue, cells or virus corresponding to Staphylococcus aureus. Soluble and structural antigens of Staphylococcus aureus whole bacterium; ATCC #27660.

### 常规说明

This antibody reacts with Staphylococcus aureus in bacterial and infected tissue samples.

The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.

If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As

### 性能

---

<b>形式</b>	Liquid
<b>存放说明</b>	Shipped at 4°C. Store at +4°C.
<b>存储溶液</b>	Preservative: 0.002% Thimerosal (merthiolate) Constituents: 1% BSA, PBS
<b>纯度</b>	Proprietary Purification
<b>Primary antibody说明</b>	This antibody reacts with Staphylococcus aureus in bacterial and infected tissue samples.
<b>克隆</b>	多克隆
<b>同种型</b>	IgG

### 应用

---

## The Abpromise guarantee

**Abpromise™** 承诺保证使用ab68952于以下的经测试应用

“应用说明”部分下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB		1/20 - 1/200.
ELISA		1/200 - 1/1000.

## 靶标

---

**相关性** S.aureus binds IgG molecules by their Fc region. In serum, the bacteria will bind IgG molecules in the wrong orientation on their surface, which hinders opsonization and phagocytosis. Mutants of S. aureus lacking protein A are more efficiently phagocytosed in vitro, and mutants in infection models have diminished virulence. Due to its affinity for the Fc region of many mammalian immunoglobulins, protein A is considered a universal reagent in biochemistry and immunology.

**细胞定位** Secreted; cell wall.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors