

HRP Anti-Trypsin Inhibitor antibody ab34606

概述

产品名称	HRP Anti-Trypsin抑制剂抗体
描述	HRP兔多克隆抗体to Trypsin抑制剂
宿主	Rabbit
偶联物	HRP
经测试应用	适用于: ELISA, Dot blot, WB
种属反应性	与反应: Soybean
免疫原	Full length protein corresponding to Soybean Trypsin Inhibitor .
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	<p>pH: 6.50</p> <p>Preservative: 0.01% Gentamicin sulphate</p> <p>Constituents: 1% BSA, 0.42% Tripotassium orthophosphate, 0.87% Sodium chloride</p>
纯度	IgG fraction
纯化说明	This antibody was purified by a multi step process which includes delipidation, salt fractionation and ion exchange chromatography followed by extensive dialysis against the buffer.
克隆	多克隆
同种型	IgG

应用

The Abpromise guarantee

Abpromise™ 承诺保证使用ab34606于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ELISA		Use at an assay dependent concentration.
Dot blot		Use at an assay dependent concentration.
WB		Use at an assay dependent concentration.

靶标

相关性

The soybean trypsin inhibitor was first crystallized by Kunitz in 1945 and is one of several such inhibitors found in soybeans. It inhibits the proteolytic action approximately of an equal weight of crystalline trypsin by combining with trypsin to form a stable compound. Chymotrypsin is only slightly inhibited by soy inhibitor.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors