

HRP Anti-Whole Human Serum antibody ab34671

概述

产品名称	HRP Anti-Whole人Serum抗体
描述	HRP兔多克隆抗体to Whole人Serum
宿主	Rabbit
偶联物	HRP
经测试应用	适用于: Electron Microscopy, ELISA, Dot blot, WB
种属反应性	与反应: Human
免疫原	Human serum proteins.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	Constituents: 1% BSA, 0.812% Sodium chloride, 0.164% Sodium phosphate
纯度	BSA is IgG and Protease free
纯化说明	IgG fraction
克隆	This product is purified by a multi-step process including delipidation, salt fractionation and ion exchange chromatography followed by extensive dialysis against buffer.
同种型	多克隆
	IgG

应用

The Abpromise guarantee

Abpromise™ 承诺保证使用ab34671于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
Electron Microscopy		
ELISA		
Dot blot		
WB		

应用说明

Dot: Use at an assay dependent dilution.
Immunoperoxidase Electron Microscopy: Use at an assay dependent dilution.
ELISA: Use at an assay dependent dilution.
WB: Use at an assay dependent dilution.
Suitable for other peroxidase-antibody based assays.

Not yet tested in other applications.
Optimal dilutions/concentrations should be determined by the end user.

靶标

相关性

The human serum proteome as a whole reflects changes in the physiological state of the human body. The onset of diseases invariably induces changes in the protein content, perhaps by producing more of one protein, less of another, or by producing a whole new protein. In all of these situations, the changes offer a direct route for diagnosing the disease.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors