

Anti-HSF1 antibody ab2923

★★★★★ [1 Abreviews](#) [13 References](#) [4 图像](#)

概述

产品名称	Anti-HSF1抗体
描述	兔多克隆抗体to HSF1
宿主	Rabbit
经测试应用	适用于: IP, WB, ICC/IF
种属反应性	与反应: Mouse, Human, Caenorhabditis elegans
免疫原	Recombinant full length protein corresponding to Human HSF1. Recombinant human HSF1 expressed in E. coli.
阳性对照	WB: HEK-293T, HeLa, K562, A431, and HepG2 whole cell lysate. C. elegans whole worm tissue. ICC/IF: HeLa and NIH3T3 cells. IP: HeLa whole cell lysate.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: 0.05% Sodium azide Constituent: 99% PBS
纯度	Whole antiserum
Primary antibody说明	All organisms respond to elevated temperatures and a variety of environmental stresses by rapid synthesis of heat shock RNAs and proteins. The regulation of heat shock gene transcription is mediated by the transcriptional activator, heat shock factor (HSF), which binds to heat shock response elements (HSEs). These HSEs are found as three repeats of a 5-nucleotide {nGAAn} module, arranged in alternating orientation and present upstream of all heat shock genes. The HSEs are highly conserved among species yet HSF purified from yeast, Drosophila and human

have different molecular weights and the proteins do not show significant immunological cross reaction. Two HSFs have been identified in human cells, HSF 1 and HSF 2, which bind to the same HSEs and have 38% sequence identity. These factors are activated by distinct stimuli, HSF 1 is responsive to classical stress signals such as heat, heavy metals and oxidative reagents, whereas HSF 2 is activated during hemin-mediated differentiation of human erythroleukemia cells. HSF 1 exists constitutively in the cytoplasm and the nucleus of unstressed cells as a monomer which lacks DNA binding activity. Through an unknown signal generated during stress, HSF 1 becomes activated to a nuclear localized, trimeric state which binds to DNA. The phosphorylation of HSF 1 is necessary for maximal transcription of heat shock genes.

克隆 多克隆
同种型 IgG

应用

The Abpromise guarantee **Abpromise™** 承诺保证使用ab2923于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
IP		Use at an assay dependent concentration. 2 µL
WB	★★★★★ (1)	1/1000 - 1/10000. Detects a band of approximately 83 kDa (predicted molecular weight: 57 kDa).
ICC/IF		1/50.

靶标

功能 DNA-binding protein that specifically binds heat shock promoter elements (HSE) and activates transcription. In higher eukaryotes, HSF is unable to bind to the HSE unless the cells are heat shocked.

序列相似性 Belongs to the HSF family.


结构域 the 9aaTAD motif is a transactivation domain present in a large number of yeast and animal transcription factors.

翻译后修饰 Phosphorylated on multiple serine residues, a subset of which are involved in stress-related regulation of transcription activation. Constitutive phosphorylation represses transcriptional activity at normal temperatures. Levels increase on specific residues heat-shock and enhance HSF1 transactivation activity. Phosphorylation on Ser-307 derepresses activation on heat-stress and in combination with Ser-303 phosphorylation appears to be involved in recovery after heat-stress. Phosphorylated on Ser-230 by CAMK2, in vitro. Cadmium also enhances phosphorylation at this site. Phosphorylation on Ser-303 is a prerequisite for HSF1 sumoylation. Phosphorylation on Ser-121 inhibits transactivation and promotes HSP90 binding. Phosphorylation on Thr-142 also mediates transcriptional activity induced by heat. Phosphorylation on Ser-326 plays an important role in heat activation of HSF1 transcriptional activity.
Sumoylated with SUMO1 and SUMO2 on heat-shock. Heat-inducible sumoylation occurs after 15 min of heat-shock, after which levels decrease and at 4 hours, levels return to control levels. Sumoylation has no effect on HSE binding nor on transcriptional activity. Phosphorylation on Ser-303 is a prerequisite for sumoylation.

细胞定位

Cytoplasm. Nucleus. Cytoplasmic during normal growth. On activation, translocates to nuclear stress granules. Colocalizes with SUMO1 in nuclear stress granules.

图片



All lanes : Anti-HSF1 antibody (ab2923) at 1/1000 dilution

Lane 1 : HEK-293T (Human epithelial cell line from embryonic kidney transformed with large T antigen) whole cell lysate

Lane 2 : HeLa (Human epithelial cell line from cervix adenocarcinoma) whole cell lysate

Lane 3 : K562 (Human chronic myelogenous leukemia cell line from bone marrow) whole cell lysate

Lane 4 : A431 (Human epidermoid carcinoma cell line) whole cell lysate

Lane 5 : HepG2 (Human liver hepatocellular carcinoma cell line) whole cell lysate

Lane 6 : COS-7 (African green monkey kidney fibroblast-like cell line) whole cell lysate

Lane 7 : NIH/3T3 (Mouse embryo fibroblast cell line) whole cell lysate

Lane 8 : NRK (Rat kidney normal tissue) whole cell lysate


Lysates/proteins at 50 µg per lane.

Secondary

All lanes : HRP-conjugated goat anti-rabbit IgG secondary antibody at 1/20000 dilution

Predicted band size: 57 kDa

Western blot analysis of Heat Shock Factor 1 (HSF1) was performed by loading samples onto a 4-20% Tris-HCl polyacrylamide gel. Proteins were transferred to a PVDF membrane and blocked with 5% BSA/TBST for at least 1 hour. The membrane was incubated with ab2923 overnight at 4°C on a rocking platform, washed in TBS-0.1% Tween 20, and probed with a secondary antibody for at least one hour. Chemiluminescent detection was performed.




Immunocytochemistry/ Immunofluorescence - Anti-HSF1 antibody (ab2923)

Immunocytochemistry/Immunofluorescence analysis of HSF1 (green) in HeLa (Human epithelial cell line from cervix adenocarcinoma) and NIH/3T3 (Mouse embryo fibroblast cell line) cells. Formalin fixed cells were permeabilized with 0.1% Triton X-100 in TBS for 10 minutes at room temperature and blocked with 1% BSA for 15 minutes at room temperature. Cells were probed with ab2923 (1:50) for at least 1 hour at room temperature, washed with PBS, and incubated with a DyLight 488-conjugated goat-anti rabbit IgG secondary antibody (1:400) for 30 minutes at room temperature. Nuclei (blue) were stained with Hoechst 33342 dye. Images were taken at 20X magnification.



Immunoprecipitation - Anti-HSF1 antibody (ab2923)

Immunoprecipitation of HSF1 was performed on HeLa (Human epithelial cell line from cervix adenocarcinoma) cells. Antigen:antibody complexes were formed by incubating 500µg whole cell lysate with 2µg of ab2923 overnight on a rocking platform at 4°C. Immune complexes were captured on 50µl Protein A/G Plus Agarose, washed extensively, and eluted with buffer. Samples were resolved on a 4-20% Tris-HCl polyacrylamide gel, transferred to a PVDF membrane, and blocked with 5% BSA/TBST for at least 1 hour. The membrane was incubated with ab2923 (1:1000) overnight rotating at 4°C, washed in TBST, and probed with detection reagent (1:1000) for at least one hour. Chemiluminescent detection was performed.



Western blot - Anti-HSF1 antibody (ab2923)

All lanes : Anti-HSF1 antibody (ab2923) at 1/5000 dilution

Lane 1 : HSF1 RNAi treated C. elegans whole worm tissue

Lane 2 : Untreated C. elegans whole worm tissue

Lysates/proteins at 10 µg per lane.

Secondary

All lanes : anti-rabbit IgG HRP-linked at 1/5000 dilution

Predicted band size: 57 kDa

Western blot analysis of HSF1 was performed by loading each sample onto an SDS-PAGE gel. Proteins were transferred to PVDF membrane and followed by blocking in TBST+5% non-fat

milk. HSF1 was detected using the primary antibody (ab2923) in TBST+2% non-fat milk overnight at 4°C on a rocking platform, followed by secondary antibody for 30 min. Detection was performed using 1-Step™ TMB-Blotting Substrate Solution.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors