

# Anti-HSV2 gB antibody [1-G-10] ab18687

## 1 References

### 概述

产品名称	Anti-HSV2 gB抗体[1-G-10]
描述	小鼠单克隆抗体[1-G-10] to HSV2 gB
宿主	Mouse
经测试应用	适用于: WB, ELISA, ICC/IF
种属反应性	与反应: Human herpesvirus 2
免疫原	Produced against HSV (II) infected cells.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	<p>pH: 7.20</p> <p>Preservative: 0.01% Sodium azide</p> <p>Constituent: PBS</p>
纯度	Protein A purified
克隆	单克隆
克隆编号	1-G-10
同种型	IgG1

### 应用

The Abpromise guarantee

**Abpromise™** 承诺保证使用ab18687于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB		Use at an assay dependent dilution.
ELISA		Use at an assay dependent dilution.
ICC/IF		Use at an assay dependent dilution.

## 靶标

**相关性** Herpes simplex type 2 (HSV-2) belongs to a family that includes HSV-1, Epstein-Barr virus (EBV) and Varicella zoster (chicken pox) virus. HSV-1 and HSV-2 are extremely difficult to distinguish from each other. These viruses have a DNA genome, an icosahedral protein coat and are encased in a lipid membrane derived from the nuclear membrane of the last host. These viruses are capable of entering a latent phase where the host shows no visible sign of infection and levels of infectious agent become very low. During the latent phase the viral DNA is integrated into the genome of the host cell.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors