

# Anti-HSV2 gG Envelope antibody [H1206] ab53470

## 1 References

### 概述

---

<b>产品名称</b>	Anti-HSV2 gG Envelope抗体[H1206]
<b>描述</b>	小鼠单克隆抗体[H1206] to HSV2 gG Envelope
<b>宿主</b>	Mouse
<b>特异性</b>	ab53470 is reactive with gG of Herpes Simplex Virus 2.
<b>经测试应用</b>	<b>适用于:</b> ICC/IF, WB
<b>种属反应性</b>	<b>与反应:</b> Human herpesvirus 2
<b>免疫原</b>	Full length protein corresponding to HSV2 gG Envelope .
<b>常规说明</b>	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

### 性能

---

<b>形式</b>	Liquid
<b>存放说明</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
<b>存储溶液</b>	pH: 7.40 Constituent: PBS
<b>纯度</b>	Protein G purified
<b>纯化说明</b>	Purified from TCS
<b>克隆</b>	单克隆
<b>克隆编号</b>	H1206
<b>同种型</b>	IgG1
<b>轻链类型</b>	kappa

## 应用

---

### The Abpromise guarantee **Abpromise™**承诺保证使用ab53470于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ICC/IF		Use a concentration of 10 µg/ml.
WB		Use a concentration of 10 µg/ml.

## 靶标

---

**相关性** HSV2 is a member of the herpes virus family, Herpesviridae. All herpes viruses are morphologically identical. They have a large double stranded DNA genome. The virion consists of an icosahedral nucleocapsid which is surrounded by a lipid bilayer envelope. Between the capsid and the envelope is an amorphous layer of proteins, termed the tegument. Glycoprotein G (gG) is one of seven envelope glycoproteins in HSV2.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors