

Anti-Human IgA antibody [3B7] ab7400

2 References

概述

产品名称	抗人IgA抗体[3B7]
描述	小鼠单克隆抗体[3B7] to人IgA
宿主	Mouse
经测试应用	适用于: ELISA, Immunodiffusion, Sandwich ELISA
种属反应性	与反应: Human
免疫原	IgA from human colostrum.
常规说明	<p>Concentration varies from lot to lot and can be provided on request.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	pH: 7.40 Preservative: 0.1% Sodium azide Constituent: PBS
纯度	Protein A purified
克隆	单克隆
克隆编号	3B7
骨髓瘤	Sp2/0
同种型	IgG1

应用

The Abpromise guarantee **Abpromise™**承诺保证使用ab7400于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ELISA		Use at an assay dependent concentration.
Immunodiffusion		Use at an assay dependent concentration.
Sandwich ELISA		Use at an assay dependent concentration. Can be paired for Sandwich ELISA with HRP Mouse monoclonal [1H9] to Human IgA (ab7383) . Can be used as Capture antibody with ab7383 and ab7402 .

靶标

相关性 Ig alpha is the major immunoglobulin class in body secretions. It may serve both to defend against local infection and to prevent access of foreign antigens to the general immunologic system.

细胞定位 Secreted

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors