


### Anti-IFIT5 antibody ab113485

#### 2 图像

#### 概述

产品名称	Anti-IFIT5抗体
描述	兔多克隆抗体to IFIT5
宿主	Rabbit
经测试应用	适用于: WB, IHC-P
种属反应性	与反应: Human 预测可用于: Mouse, Rat 
免疫原	Recombinant fragment, corresponding to amino acids 79-294 of Human IFIT5 (BC025786).
阳性对照	Human fetal brain and fetal colon lysates; Human fetal brain tissue.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

#### 性能

形式	Lyophilized:Reconstitute in 200ul sterile H2O.
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
存储溶液	pH: 7.20 Preservative: 0.02% Sodium azide Constituent: 1% BSA
纯度	Immunogen affinity purified
纯化说明	ab113485 is purified by a protein affinity column.
克隆	多克隆
同种型	IgG

#### 应用

## The Abpromise guarantee

**Abpromise™** 承诺保证使用ab113485于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。


应用	Ab评论	说明
WB		1/1000 - 1/2000. Detects a band of approximately 50 kDa (predicted molecular weight: 56 kDa).
IHC-P		1/100 - 1/500.

## 靶标

### 序列相似性


Belongs to the IFIT family.  
Contains 8 TPR repeats.

## 图片



ab113485 at 1/100 dilution staining IFIT5 in Formalin-fixed, Paraffin-embedded Human fetal brain tissue by Immunohistochemistry.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-IFIT5 antibody (ab113485)



**All lanes :** Anti-IFIT5 antibody (ab113485) at 1/1000 dilution

**Lane 1 :** Human fetal brain lysate

**Lane 2 :** Human fetal colon lysate

**Predicted band size:** 56 kDa

Western blot - Anti-IFIT5 antibody (ab113485)

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet

- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors