

Anti-IgA2 antibody ab88250

4 图像

概述

产品名称	Anti-IgA2抗体
描述	小鼠多克隆抗体to IgA2
宿主	Mouse
经测试应用	适用于: WB, ICC/IF
种属反应性	与反应: Human
免疫原	Recombinant full length protein within Human IgA2. The exact immunogen sequence used to generate this antibody is proprietary information. If additional detail on the immunogen is needed to determine the suitability of the antibody for your needs, please contact our Scientific Support team to discuss your requirements.
阳性对照	WB: Human colon tissue lysate; Hela nuclear extract; IgA2 transfected 293T cell lysate ICC/IF: HeLa cell
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	pH: 7.40 Constituent: 100% PBS
纯度	Protein A purified
克隆	多克隆
同种型	IgG

应用

The Abpromise guarantee Abpromise™ 承诺保证使用ab88250于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。


应用	Ab评论	说明
WB		Use a concentration of 1 µg/ml. Predicted molecular weight: 53 kDa.
ICC/IF		Use a concentration of 10 µg/ml.

靶标

相关性 Ig alpha is the major immunoglobulin class in body secretions. It may serve both to defend against local infection and to prevent access of foreign antigens to the general immunologic system.

细胞定位 Secreted

图片



Anti-IgA2 antibody (ab88250) at 1 µg/ml + human colon tissue lysate at 50 µg

Predicted band size: 53 kDa

Observed band size: 53 kDa




Anti-IgA2 antibody (ab88250) at 1 µg/ml + Hela nuclear extract at 50 µg

Predicted band size: 53 kDa

Observed band size: 35 kDa

Additional bands at: 48 kDa. We are unsure as to the identity of these extra bands.



Western blot - Anti-IgA2 antibody (ab88250)

All lanes : Anti-IgA2 antibody (ab88250) at 1 µg/ml


Lane 1 : IgA2 transfected 293T cell lysate

Lane 2 : Non transfected 293T cell lysate

Lysates/proteins at 50 µg per lane.

Predicted band size: 53 kDa

Observed band size: 53 kDa



Immunocytochemistry/ Immunofluorescence - Anti-IgA2 antibody (ab88250)

ab88250, at 10 µg/ml, staining IgA2 in HeLa cell by Immunofluorescence.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors

