

Anti-IL-13 antibody [A130D 12G5 1E4] ab52538

概述

产品名称	Anti-IL-13抗体[A130D 12G5 1E4]
描述	小鼠单克隆抗体[A130D 12G5 1E4] to IL-13
宿主	Mouse
经测试应用	适用于: ELISA, Dot blot, Flow Cyt, WB
种属反应性	与反应: Rat
免疫原	Recombinant full length protein corresponding to Rat IL-13.
阳性对照	FACS: intracellular IL-13 can be detected in PMA/ionophore-stimulated rat splenocytes
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term.
存储溶液	pH: 7.20 Constituent: PBS
纯度	Protein G purified
纯化说明	>95% pure as determined by SDS-PAGE analysis. 0.22 micron sterile filtered.
克隆	单克隆
克隆编号	A130D 12G5 1E4
同种型	IgG1

应用

The Abpromise guarantee

Abpromise™ 承诺保证使用ab52538于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ELISA		Use at an assay dependent concentration.
Dot blot		Use at an assay dependent concentration.
Flow Cyt		Use at an assay dependent concentration. ab170190 - Mouse monoclonal IgG1, is suitable for use as an isotype control with this antibody.
WB		Use at an assay dependent concentration.

靶标

功能	Cytokine (PubMed:8096327, PubMed:8097324). Inhibits inflammatory cytokine production (PubMed:8096327). Synergizes with IL2 in regulating interferon-gamma synthesis (PubMed:8096327). May be critical in regulating inflammatory and immune responses (PubMed:8096327, PubMed:8097324). Positively regulates IL31RA expression in macrophages.
疾病相关	Allergic rhinitis
序列相似性	Belongs to the IL-4/IL-13 family.
细胞定位	Secreted.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors