

# Anti-IL-17B antibody ab106272

## 2 图像

### 概述

产品名称	Anti-IL-17B抗体
描述	兔多克隆抗体to IL-17B
宿主	Rabbit
经测试应用	适用于: WB, IHC-P
种属反应性	与反应: Human 预测可用于: Mouse, Rat 
免疫原	Recombinant full length protein corresponding to Human IL-17B aa 21-180. Database link: <a href="#">BC113925</a>
阳性对照	HeLa cell lysate; Human fetal small intestine tissue.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

### 性能

形式	Lyophilized:Reconstitute in 200ul sterile water
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	Preservative: 0.02% Sodium azide
纯度	Protein A purified
克隆	多克隆
同种型	IgG

### 应用

**The Abpromise guarantee** **Abpromise™**承诺保证使用ab106272于以下的经测试应用


“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB		1/500 - 1/1000. Detects a band of approximately 20 kDa (predicted molecular weight: 20 kDa).
IHC-P		1/100 - 1/500.

## 靶标


功能	Stimulates the release of tumor necrosis factor alpha and IL-1-beta from the monocytic cell line THP-1.
组织特异性	Expressed in adult pancreas, small intestine, stomach, spinal cord and testis. Less pronounced expression in prostate, colon mucosal lining, and ovary.
序列相似性	Belongs to the IL-17 family.
细胞定位	Secreted.

## 图片



Anti-IL-17B antibody (ab106272) at 1/500 dilution + HeLa cell lysate

**Predicted band size:** 20 kDa



ab106272 at 1/100 dilution staining IL17B in Human fetal small intestine by Immunohistochemistry, Formalin-fixed, Paraffin-embedded tissue.

## **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors