

Anti-IL-20R2 antibody [20RNTC] ab95824

[5 References](#) [2 图像](#)

概述

产品名称	Anti-IL-20R2抗体[20RNTC]
描述	大鼠单克隆抗体[20RNTC] to IL-20R2
宿主	Rat
特异性	ab95824 recognizes Mouse and Human IL-20R2 in either of its two heterodimeric complexes formed by pairing with IL-20RA or IL-22RA1
经测试应用	适用于: Flow Cyt, ICC/IF
种属反应性	与反应: Mouse, Human
免疫原	Recombinant fragment corresponding to Mouse IL-20R2 (extracellular).
阳性对照	Fixed and permeabilized A375 cells
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	pH: 7.20 Preservative: 0.09% Sodium azide Constituent: 1% PBS
纯度	Protein G purified
克隆	单克隆
克隆编号	20RNTC
同种型	IgG2a
轻链类型	kappa

应用

The Abpromise guarantee **Abpromise™**承诺保证使用ab95824于以下的经测试应用


“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
Flow Cyt		Use a concentration of 1 - 10 µg/ml. ab18450 - Rat monoclonal IgG2a, is suitable for use as an isotype control with this antibody.
ICC/IF		Use a concentration of 1 - 10 µg/ml.


靶标

功能	The IL20RA/IL20RB dimer is a receptor for IL19, IL20 and IL24. The IL22RA1/IL20RB dimer is a receptor for IL20 and IL24.
组织特异性	Widely expressed with highest levels in skin and testis. Highly expressed in psoriatic skin.
序列相似性	Belongs to the type II cytokine receptor family. Contains 2 fibronectin type-III domains.
细胞定位	Membrane.

图片



Immunofluorescent microscopy of fixed and permeabilized A375 cells using 10 µg/ml of ab95824 (right) or Rat IgG2a - Isotype Control Purified (left) followed by Anti-Rat IgG Biotin and then Streptavidin FITC.



Mouse IL-20R2-GFP 293 transfectants stained with 0.25ug ab95824 + anti-rat IgG PE

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors