

# Anti-ILT-7 antibody [17G10.2] ab234070

[1 References](#) [1 图像](#)

### 概述

<b>产品名称</b>	Anti-ILT-7抗体[17G10.2]
<b>描述</b>	小鼠单克隆抗体[17G10.2] to ILT-7
<b>宿主</b>	Mouse
<b>经测试应用</b>	<b>适用于:</b> Flow Cyt
<b>种属反应性</b>	<b>与反应:</b> Human
<b>免疫原</b>	The details of the immunogen for this antibody are not available.
<b>表位</b>	extracellular
<b>阳性对照</b>	Flow Cytometry: Human buffy coat cells.
<b>常规说明</b>	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

### 性能

<b>形式</b>	Liquid
<b>存放说明</b>	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
<b>存储溶液</b>	pH: 7.40 Preservative: 0.0975% Sodium azide Constituent: PBS
<b>纯度</b>	Protein A purified
<b>纯化说明</b>	Purity >95% by SDS-PAGE.
<b>克隆</b>	单克隆
<b>克隆编号</b>	17G10.2
<b>同种型</b>	IgG1
<b>轻链类型</b>	kappa

## 应用

**The Abpromise guarantee**      **Abpromise™**承诺保证使用ab234070于以下的经测试应用


“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
Flow Cyt		Use a concentration of 2 - 5 µg/ml.

## 靶标

<b>功能</b>	May act as receptor for class I MHC antigens. Ligand binding leads to the activation of eosinophils and to the release of RNASE2, IL4 and leukotriene C4.
<b>组织特异性</b>	Detected in eosinophils, neutrophils and monocytes.
<b>序列相似性</b>	Contains 4 Ig-like C2-type (immunoglobulin-like) domains.
<b>细胞定位</b>	Membrane.

## 图片



Surface staining of ILT-7 in human buffy coat cells with ab234070 followed by a Goat anti-mouse-APC conjugated secondary.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors