

Anti-Integrin alpha 4+beta 7 antibody [RM0059-2G18] ab73261

1 图像

概述	
产品名称	Anti-Integrin alpha 4+beta 7抗体[RM0059-2G18]
描述	大鼠单克隆抗体[RM0059-2G18] to Integrin alpha 4+beta 7
宿主	Rat
经测试应用	适用于: Flow Cyt, IP, IHC-P
种属反应性	与反应: Mouse
免疫原	Cell lysates of mouse lymphocytes
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>
性能	
形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Constituent: PBS
纯度	Protein A purified
克隆	单克隆
克隆编号	RM0059-2G18
同种型	IgG2a
应用	

The Abpromise guarantee

Abpromise™承诺保证使用ab73261于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。


应用	Ab评论	说明
Flow Cyt		Use at an assay dependent concentration.
IP		Use at an assay dependent concentration.
IHC-P		1/100 - 1/200.

靶标

相关性 Integrin alpha 4 beta 7, also known as, the lymphocyte Peyer's patch adhesion molecule (LPAM1) is a member of the integrin family of cell surface receptors. It is expressed primarily on mucosal lymphocytes, but is also present on NK cells and eosinophils. The alpha 4 beta 7 heterodimer mediates the binding of lymphocytes to its ligand, mucosal vascular addressin (MAdCAM1) on the high endothelial venules, thereby directing the homing of lymphocytes into Peyer's patches and intestinal lamina propria.

细胞定位 Membrane; Single pass type I membrane protein.

图片



ab73261 at 1/200 dilution, staining Integrin alpha 4+beta 7 in LPS exposed mouse spleen tissue section by Immunohistochemistry (Formalin/PFA fixed paraffin-embedded sections).

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Integrin alpha 4+beta 7 antibody [RM0059-2G18] (ab73261)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish

- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors