

Anti-Interferon alpha 1 antibody [ST29] ab20200

[2 References](#) [1 图像](#)

概述

产品名称	Anti-Interferon alpha 1 抗体[ST29]
描述	小鼠单克隆抗体[ST29] to Interferon alpha 1
宿主	Mouse
经测试应用	适用于: Sandwich ELISA
种属反应性	与反应: Recombinant fragment
免疫原	Human interferon (HuIFN) produced by Namalwa cells (HuIFN-aN).
常规说明	Myeloma: NS2 The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing. If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: 0.02% Sodium azide Constituent: 99.98% PBS
纯度	Protein A/G purified
克隆	单克隆
克隆编号	ST29
同种型	IgG1
轻链类型	unknown

应用

The Abpromise guarantee

Abpromise™承诺保证使用ab20200于以下的经测试应用


“应用说明”部分下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
Sandwich ELISA		Use a concentration of 5 µg/ml.

靶标

功能	Produced by macrophages, IFN-alpha have antiviral activities. Interferon stimulates the production of two enzymes: a protein kinase and an oligoadenylate synthetase.
序列相似性	Belongs to the alpha/beta interferon family.
细胞定位	Secreted.

图片



Standard Curve for Interferon alpha (Analyte: **ab48750**) dilution range 1pg/ml to 1ug/ml using Capture Antibody Mouse monoclonal [ST29] to Interferon alpha (ab20200) at 5ug/ml and Detector Antibody Rabbit polyclonal to Interferon alpha (Biotin) (**ab84257**) at 0.5ug/ml.

Sandwich ELISA - Anti-Interferon alpha 1 antibody
[ST29] (ab20200)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors