

Anti-Interferon alpha 1 antibody [UMCh3-IFN-a] ab24780

概述

产品名称	Anti-Interferon alpha 1 抗体[UMCh3-IFN-a]
描述	小鼠单克隆抗体[UMCh3-IFN-a] to Interferon alpha 1
宿主	Mouse
特异性	ab24780 recognises interferon alpha 1, does not cross react with IFN beta or gamma and does not neutralize chicken IFN alpha.
经测试应用	适用于: ELISA, WB
种属反应性	与反应: Chicken
免疫原	Recombinant full length protein corresponding to Chicken Interferon alpha 1.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term.
存储溶液	Preservative: 0.09% Sodium azide
纯度	Protein G purified
克隆	单克隆
克隆编号	UMCh3-IFN-a
同种型	IgG1

应用

The Abpromise guarantee

Abpromise™ 承诺保证使用ab24780于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ELISA		Use at an assay dependent concentration.
WB		Use at an assay dependent concentration.

靶标

功能	Produced by macrophages, IFN-alpha have antiviral activities. Interferon stimulates the production of two enzymes: a protein kinase and an oligoadenylate synthetase.
序列相似性	Belongs to the alpha/beta interferon family.
细胞定位	Secreted.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors