

Anti-Interferon gamma antibody [EP1109Y] ab52746

RabMAb


1 References [2 图像](#)

概述

产品名称	Anti-Interferon gamma抗体[EP1109Y]
描述	兔单克隆抗体[EP1109Y] to Interferon gamma
宿主	Rabbit
经测试应用	适用于: WB 不适用于: ICC/IF or IHC-P
种属反应性	与反应: Human, Recombinant fragment
免疫原	Synthetic peptide. This information is proprietary to Abcam and/or its suppliers.
阳性对照	WB: Human Interferon gamma recombinant protein.
常规说明	<p>Our RabMAb[®] technology is a patented hybridoma-based technology for making rabbit monoclonal antibodies. For details on our patents, please refer to RabMAb[®] patents.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p> <p>Mouse, Rat: We have preliminary internal testing data to indicate this antibody may not react with these species. Please contact us for more information.</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C. Avoid freeze / thaw cycle.
存储溶液	pH: 7.20 Preservative: 0.01% Sodium azide Constituents: 59% PBS, 40% Glycerol (glycerin, glycerine), 0.15% BSA
纯度	Protein A purified



Western blot - Anti-Interferon gamma antibody
[EP1109Y] (ab52746)

Anti-Interferon gamma antibody [EP1109Y] (ab52746) at 1/1000 dilution + Human Interferon gamma recombinant protein at 0.015 µg


Secondary

Goat Anti-Rabbit IgG H&L (HRP) ([ab97051](#)) at 1/20000 dilution

Predicted band size: 19 kDa

Blocking buffer: 5% NFDM/TBST

Diluting buffer: 5% NFDM/TBST



Western blot - Anti-Interferon gamma antibody
[EP1109Y] (ab52746)

Anti-Interferon gamma antibody [EP1109Y] (ab52746) at 1/2000 dilution + recombinant protein at 0.01 µg

Secondary

Goat anti-rabbit HRP labeled at 1/2000 dilution

Predicted band size: 19 kDa

Observed band size: 17 kDa

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet

- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors