

Anti-IRF2BP2 antibody ab114997

1 图像

概述

产品名称	Anti-IRF2BP2抗体
描述	兔多克隆抗体to IRF2BP2
宿主	Rabbit
经测试应用	适用于: IP
种属反应性	与反应: Human 预测可用于: Chimpanzee, Rhesus monkey 
免疫原	Synthetic peptide, corresponding to a region within amino acids 225-275 of Human IRF2BP2 (NP_892017.2).
阳性对照	Whole cell lysate from Hela cells
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
存储溶液	pH: 7 Preservative: 0.09% Sodium azide Constituent: 99% Tris citrate/phosphate
纯度	pH 7 to 8 Immunogen affinity purified
纯化说明	ab114997 was affinity purified using an epitope specific to IRF2BP2 immobilized on solid support.
克隆	多克隆

应用

The Abpromise guarantee **Abpromise™**承诺保证使用ab114997于以下的经测试应用


“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
IP		Use at 2-10 µg/mg of lysate.

靶标

功能	Acts as a transcriptional repressor. Acts as a transcriptional corepressor in a IRF2-dependent manner. This repression is not mediated at least in part by histone deacetylase activities.
序列相似性	Belongs to the IRF2BP family.
结构域	The C-terminal RING-type zinc finger domain is necessary for interaction with BAP1.
细胞定位	Nucleus.

图片



Immunoprecipitation: 1 mg of HeLa cell lysate immunoprecipitated using ab114997 at 6 µg/mg lysate; 20% of immunoprecipitate was loaded/lane.

Western Blot of immunoprecipitate: Anti-IRF2BP2 antibody at 1 µg/ml.

Developed using the ECL technique

Exposure time: 1 second

Predicted band size: 61 kDa

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise,

please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors