

# Anti-Kanamycin Receptor antibody ab23971

### 概述

|       |   |
|-------|---|
| 产品名称  | Anti-Kanamycin Receptor抗体   |
| 描述    | 兔多克隆抗体to Kanamycin  |
| 宿主    | Rabbit  |
| 特异性   | Antiserum reacts specifically with kanamycin.   |
| 经测试应用 | 适用于: WB, IHC, ELISA   |
| 种属反应性 | 与反应: Species independent  |
| 免疫原   | Recombinant full length protein corresponding to Kanamycin.   |
| 常规说明  | <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p> |

### 性能

|      |   |
|------|---|
| 形式   | Liquid  |
| 存放说明 | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. |
| 存储溶液 | Constituent: Whole serum  |
| 纯度   | Whole antiserum   |
| 克隆   | 多克隆   |
| 同种型  | IgG   |

### 应用

**The Abpromise guarantee**      **Abpromise™** 承诺保证使用ab23971于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

| 应用    | Ab评论 | 说明                 |
|-------|------|--------------------|
| WB    |      | 1/200 - 1/1000.    |
| IHC   |      | 1/100 - 1/500.     |
| ELISA |      | 1/15000 - 1/50000. |

## 靶标

**相关性** Kanamycin is a water-soluble aminoglycoside antibiotic that is derived from the bacterium *Streptomyces kanamyceticus*. The drug has a very low adsorption rate in the intestine and is used only for topical effect.

**细胞定位** Secreted

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors