

Anti-KAP1 antibody [20C1] ab22553

★★★★★ [14 Abreviews](#) [78 References](#) [4 图像](#)

概述

产品名称	Anti-KAP1抗体[20C1]
描述	小鼠单克隆抗体[20C1] to KAP1
宿主	Mouse
经测试应用	适用于: WB, ICC, IHC-P
种属反应性	与反应: Human
免疫原	Recombinant fragment, corresponding to amino acids 60-383 of Human KAP1.
阳性对照	IHC-P: Human spleen tissue. WB: HEK-293 (untreated and treated with hydrogen peroxide at 100 uM for 10 min), A549, HEL 92.1.7, PC-3, Caco-2 and BJ fibroblast cell lysates. ICC: HeLa cells.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: 0.05% Sodium azide Constituents: PBS, 0.1% BSA
纯度	Protein G purified
克隆	单克隆
克隆编号	20C1
同种型	IgG1

应用

The Abpromise guarantee
Abpromise™承诺保证使用ab22553于以下的经测试应用
“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB	★★★★★ (7)	1/500. Predicted molecular weight: 88 kDa.
ICC		Use a concentration of 2 µg/ml.
IHC-P		Use a concentration of 5 µg/ml.

靶标	
功能	<p>Nuclear corepressor for KRAB domain-containing zinc finger proteins (KRAB-ZFPs). Mediates gene silencing by recruiting CHD3, a subunit of the nucleosome remodeling and deacetylation (NuRD) complex, and SETDB1 (which specifically methylates histone H3 at 'Lys-9' (H3K9me)) to the promoter regions of KRAB target genes. Enhances transcriptional repression by coordinating the increase in H3K9me, the decrease in histone H3 'Lys-9 and 'Lys-14' acetylation (H3K9ac and H3K14ac, respectively) and the disposition of HP1 proteins to silence gene expression. Recruitment of SETDB1 induces heterochromatinization. May play a role as a coactivator for CEBPB and NR3C1 in the transcriptional activation of ORM1. Also corepressor for ERBB4. Inhibits E2F1 activity by stimulating E2F1-HDAC1 complex formation and inhibiting E2F1 acetylation. May serve as a partial backup to prevent E2F1-mediated apoptosis in the absence of RB1. Important regulator of CDKN1A/p21(CIP1). Has E3 SUMO-protein ligase activity toward itself via its PHD-type zinc finger.</p>
组织特异性	<p>Expressed in all tissues tested including spleen, thymus, prostate, testis, ovary, small intestine, colon and peripheral blood leukocytes.</p>
通路	<p>Protein modification; protein sumoylation.</p>
序列相似性	<p>Belongs to the TRIM/RBCC family.</p> <p>Contains 2 B box-type zinc fingers.</p> <p>Contains 1 bromo domain.</p> <p>Contains 1 PHD-type zinc finger.</p> <p>Contains 1 RING-type zinc finger.</p>
结构域	<p>The HP1 box is both necessary and sufficient for HP1 binding.</p> <p>The PHD-type zinc finger enhances CEBPB transcriptional activity. The PHD-type zinc finger, the HP1 box and the bromo domain, function together to assemble the machinery required for repression of KRAB domain-containing proteins. Acts as an intramolecular SUMO E3 ligase for autosumoylation of bromodomain.</p> <p>The RING-finger-B Box-coiled-coil/tripartite motif (RBCC/TRIM motif) is required for interaction with the KRAB domain of KRAB-zinc finger proteins. Binds four zinc ions per molecule. The RING finger and the N-terminal of the leucine zipper alpha helical coiled-coil region of RBCC are required for oligomerization.</p> <p>Contains one Pro-Xaa-Val-Xaa-Leu (PxVxL) motif, which is required for interaction with chromoshadow domains. This motif requires additional residues -7, -6, +4 and +5 of the central Val which contact the chromoshadow domain.</p>
翻译后修饰	<p>Phosphorylated upon DNA damage, probably by ATM or ATR. ATM-induced phosphorylation on Ser-824 represses sumoylation leading to the de-repression of expression of a subset of genes involved in cell cycle control and apoptosis in response to genotoxic stress. Dephosphorylation by the phosphatases, PPP1CA and PP1CB forms, allows sumoylation and expression of TRIM28</p>


target genes.

Sumoylation/desumoylation events regulate TRIM28-mediated transcriptional repression. Sumoylation is required for interaction with CHD3 and SETDB1 and the corepressor activity. Represses and is repressed by Ser-824 phosphorylation. Enhances the TRIM28 corepressor activity, inhibiting transcriptional activity of a number of genes including GADD45A and CDKN1A/p21. Lys-554, Lys-779 and Lys-804 are the major sites of sumoylation. In response to Dox-induced DNA damage, enhanced phosphorylation on Ser-824 prevents sumoylation and allows de-repression of CDKN1A/p21.

细胞定位

Nucleus. Associated with centromeric heterochromatin during cell differentiation through CBX1.

图片



All lanes : Anti-KAP1 antibody [20C1] (ab22553) at 1 µg/ml

Lanes 1-2 : HEK-293 (Human epithelial cell line from embryonic kidney) whole cell lysate

Lane 3 : A549 (Human lung carcinoma cell line) whole cell lysate

Lane 4 : HEL 92.1.7 whole cell lysate

Lane 5 : PC-3 (Human prostate adenocarcinoma cell line) whole cell lysate

Lane 6 : Caco-2 (Human colorectal adenocarcinoma cell line) whole cell lysate


Lysates/proteins at 30 µg per lane.

Predicted band size: 88 kDa

Observed band size: 110 kDa


Western blot analysis was performed on nuclear enriched extracts (30 µg lysate) of cells treated with Hydrogen peroxide (100 µM Hydrogen peroxide for 10 min) or untreated. The blots were probed with ab22553 at 1 µg/mL and detected by chemiluminescence using Goat anti-Mouse IgG (H+L) Superclonal™ Secondary Antibody, HRP conjugate. A 110 kDa band corresponding to KAP1 was observed across the cell lines tested. Known quantity of protein samples were electrophoresed using Novex® NuPAGE® 4-12 % Bis-Tris gel, XCell SureLock™ Electrophoresis System and Novex® Sharp Pre-Stained Protein Standard. Resolved proteins were then transferred onto a nitrocellulose membrane with iBlot® 2 Dry Blotting System. The membrane was probed with the relevant primary and secondary Antibody following blocking with 5 % skimmed milk. Chemiluminescent detection was performed using

Pierce™ ECL Western Blotting Substrate.



Immunocytochemistry - Anti-KAP1 antibody [20C1]
(ab22553)


Immunofluorescence analysis of TIF1 beta was performed using 70% confluent log phase HeLa (Human cervix adenocarcinoma epithelial cell) cells. The cells were fixed with 4% paraformaldehyde for 10 minutes, permeabilized with 0.1% Triton™ X-100 for 10 minutes, and blocked with 1% BSA for 1 hour at room temperature. The cells were labeled with ab22553 at 2 µg/mL in 0.1% BSA and incubated for 3 hours at room temperature and then labeled with Goat anti-Mouse IgG (H+L) Superclonal™ Secondary Antibody, Alexa Fluor® 488 conjugate at a dilution of 1:2000 for 45 minutes at room temperature (Panel a: green). Nuclei (Panel b: blue) were stained with SlowFade® Gold Antifade Mountant with DAPI. F-actin (Panel c: red) was stained with Rhodamine Phalloidin. Panel d represents the merged image showing nuclear localization. Panel e shows the no primary antibody control. The images were captured at 60X magnification.



Immunohistochemistry (Formalin/PFA-fixed paraffin-
embedded sections) - Anti-KAP1 antibody [20C1] -
ChIP Grade (ab22553)

IHC image of KAP1 staining in Human normal spleen formalin fixed paraffin embedded tissue section, performed on a Leica Bond™ system using the standard protocol F. The section was pre-treated using heat mediated antigen retrieval with sodium citrate buffer (pH6, epitope retrieval solution 1) for 20 mins. The section was then incubated with ab22553, 5µg/ml, for 15 mins at room temperature and detected using an HRP conjugated compact polymer system. DAB was used as the chromogen. The section was then counterstained with haematoxylin and mounted with DPX.

For other IHC staining systems (automated and non-automated) customers should optimize variable parameters such as antigen retrieval conditions, primary antibody concentration and antibody incubation times.



Anti-KAP1 antibody [20C1] (ab22553) at 1/500 dilution + Whole cell human BJ fibroblast protein lysate at 10 µg

Predicted band size: 88 kDa

Western blot analysis of TRIM28 was performed by loading 10 µg of whole cell human BJ fibroblast protein lysate and run on a 4-12% BTE gel. Proteins were transferred to PVDF membrane. Membrane was blocked in 5% non-fat milk in TBST. TRIM28 was detected at approximately 90 kDa using ab22553 at a dilution of 1:500 (2 µg/mL) in 5% milk, followed by a 1:10,000 dilution of anti-mouse HRP. Chemiluminescent detection was performed using SuperSignal West Pico PLUS substrate.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors