# abcam

### Product datasheet

## Anti-KCTD20 antibody ab122094

1 References 5 图像

概述

产**品名称** Anti-KCTD20抗体

描述 兔多克隆抗体to KCTD20

**宿主** Rabbit

经测试应用 适用于: ICC/IF, IHC-P

种属反应性 与反应: Human

免疫原 antigen sequence: LAEDIKGSCF QSGNKRNHEP FIAPERFGNS SVGFGSNSHS

QAPEKVTLLV DGTRFVVNPQ IFTAHPDTML, corresponding to amino acids 74-143 of Human

KCTD20.

Run BLAST with EXPASY Run BLAST with S NCBI

阳性对照 ICC/IF: U-251 MG cell lysates. IHC: Human prostate, kidney, cervix and placenta tissue.

常规说明

The Life Science industry has been in the grips of a reproducibility crisis for a number of years.

Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets

your needs before purchasing.

If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be

found below, along with publications, customer reviews and Q&As

性能

形式 Liquid

存放说明 Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

**存储溶液** pH: 7.20

Preservative: 0.02% Sodium azide

Constituents: 59% PBS, 40% Glycerol (glycerin, glycerine)

纯**度** Immunogen affinity purified

**克隆** 多克隆

**同种型** lgG

应用

1

#### The Abpromise guarantee

#### Abpromise™承诺保证使用ab122094于以下的经测试应用

"应用说明"部分下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。


应用	Ab评论	说明
ICC/IF		Use a concentration of 0.25 - 2 µg/ml.
IHC-P		1/50 - 1/200. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.

#### 靶标

序列相似性


Contains 1 BTB (POZ) domain.

#### 图片



Immunocytochemistry/ Immunofluorescence - Anti-KCTD20 antibody (ab122094)

Immunofluorescent staining of Human cell line U-251MG shows positivity in endoplasmic reticulum. Recommended concentration of ab122094 1-4  $\mu$ g/ml. Cells treated with PFA/Triton X-100.




Immunohistochemistry (Formalin/PFA-fixed paraffinembedded sections) - Anti-KCTD20 antibody (ab122094)

ab122094, at 1/35 dilution, staining KCTD20 in paraffin-embedded Human kidney tissue by Immunohistochemistry.




Immunohistochemistry (Formalin/PFA-fixed paraffinembedded sections) - Anti-KCTD20 antibody (ab122094)

ab122094, at 1/35 dilution, staining KCTD20 in paraffin-embedded Human cervix tissue by Immunohistochemistry.



Immunohistochemistry (Formalin/PFA-fixed paraffinembedded sections) - Anti-KCTD20 antibody (ab122094)

ab122094, at 1/35 dilution, staining KCTD20 in paraffin-embedded Human placenta tissue by Immunohistochemistry.



Immunohistochemistry (Formalin/PFA-fixed paraffinembedded sections) - Anti-KCTD20 antibody (ab122094)

ab122094, at 1/35 dilution, staining KCTD20 in paraffin-embedded Human prostate tissue by Immunohistochemistry.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

#### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise,

please visit <a href="https://www.abcam.cn/abpromise">https://www.abcam.cn/abpromise</a> or contact our technical team.

#### Terms and conditions

• Guarantee only valid for products bought direct from Abcam or one of our authorized distributors