

Anti-KDM1/LSD1 antibody [EPR6825] - Nuclear Marker and ChIP Grade ab129195

敲除验证
重组
RabMAb

★★★★☆
[5 Abreviews](#)
[26 References](#)
[19 图像](#)

概述

产品名称	Anti-KDM1/LSD1抗体[EPR6825] -核Marker and ChIP Grade
描述	兔单克隆抗体[EPR6825] to KDM1/LSD1 -核Marker and ChIP Grade
宿主	Rabbit
经测试应用	适用于: Flow Cyt (Intra), ChIC/CUT&RUN-seq, WB, IP, IHC-P, ICC/IF, ChIP
种属反应性	与反应: Mouse, Rat, Human
免疫原	Synthetic peptide. This information is proprietary to Abcam and/or its suppliers.
阳性对照	WB: HAP1, 293T, HEK293, HeLa, Jurkat, PC3, C6, Raw 264.7, PC-12, and NIH 3T3 cell lysates. IHC-P: Human testis, rat kidney, and mouse colon tissues. ICC/IF: HAP1 and HeLa cells. Flow Cyt (intra): HeLa cells. IP: Jurkat cell lysate. ChIP: HCT 116 (Human colorectal carcinoma epithelial cell). ChIC/CUT&RUN-Seq: HeLa cells.
常规说明	<p>This product is a recombinant monoclonal antibody, which offers several advantages including:</p> <ul style="list-style-type: none"> - High batch-to-batch consistency and reproducibility - Improved sensitivity and specificity - Long-term security of supply - Animal-free production <p>For more information see here.</p> <p>Our RabMAb[®] technology is a patented hybridoma-based technology for making rabbit monoclonal antibodies. For details on our patents, please refer to RabMAb[®] patents.</p>

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C. Avoid freeze / thaw cycle.
存储溶液	<p>pH: 7.20</p> <p>Preservative: 0.01% Sodium azide</p> <p>Constituents: PBS, 40% Glycerol, 0.05% BSA</p>
纯度	Protein A purified

克隆	单克隆
克隆编号	EPR6825
同种型	IgG

应用

The Abpromise guarantee

Abpromise™承诺保证使用ab129195于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。


应用	Ab评论	说明
Flow Cyt (Intra)		1/20. ab172730 - Rabbit monoclonal IgG, is suitable for use as an isotype control with this antibody.
ChIC/CUT&RUN-seq		Use at an assay dependent concentration. 5µg
WB	★★★★★ (2)	1/10000 - 1/50000. Detects a band of approximately 110 kDa (predicted molecular weight: 92 kDa).
IP		1/10 - 1/100.
IHC-P		1/50 - 1/250. Perform heat mediated antigen retrieval before commencing with IHC staining protocol.
ICC/IF	★★★★★ (2)	Use a concentration of 1 µg/ml.
ChIP	★★★★★ (1)	Use 5 µg for 25 µg of chromatin.

靶标

功能	Histone demethylase that demethylates both 'Lys-4' (H3K4me) and 'Lys-9' (H3K9me) of histone H3, thereby acting as a coactivator or a corepressor, depending on the context. Acts by oxidizing the substrate by FAD to generate the corresponding imine that is subsequently hydrolyzed. Acts as a corepressor by mediating demethylation of H3K4me, a specific tag for epigenetic transcriptional activation. Demethylates both mono- (H3K4me1) and di-methylated (H3K4me2) H3K4me. May play a role in the repression of neuronal genes. Alone, it is unable to demethylate H3K4me on nucleosomes and requires the presence of RCOR1/CoREST to achieve such activity. Also acts as a coactivator of androgen receptor (ANDR)-dependent transcription, by being recruited to ANDR target genes and mediating demethylation of H3K9me, a specific tag for epigenetic transcriptional repression. The presence of PRKCB in ANDR-containing complexes, which mediates phosphorylation of 'Thr-6' of histone H3 (H3T6ph), a specific tag that prevents demethylation H3K4me, prevents H3K4me demethylase activity of KDM1A. Demethylates di-methylated 'Lys-370' of p53/TP53 which prevents interaction of p53/TP53 with TP53BP1 and represses p53/TP53-mediated transcriptional activation. Demethylates and stabilizes the DNA methylase DNMT1. Required for gastrulation during embryogenesis.
组织特异性	Ubiquitously expressed.
序列相似性	Belongs to the flavin monoamine oxidase family.

	Contains 1 SWIRM domain.
结构域	The SWIRM domain may act as an anchor site for a histone tail.
细胞定位	Nucleus.


图片



ChIC/CUT&RUN was performed using a pAG-MNase at a final concentration of 700 ng/mL. 2.5×10^5 of Human wild-type HeLa cell line (**ab255928**) or KDM1A (LSD1) knockout HeLa cell line (**ab265790**) were used along with 5µg of ab129195 [EPR6825]. Assay Quality Control was conducted using 5µg Anti-CTCF(**ab188408**) on the same cell lines. The resulting DNA was sequenced on the Illumina NovaSeq 6000 to a depth of 10 million reads. The negative IgG control **ab172730** is also shown.

Additional screenshots of mapped reads can be downloaded [here](#).

The University of Geneva owns patents relevant to ChIC (Chromatin Immuno-Cleavage) methods.



All lanes : Anti-KDM1/LSD1 antibody [EPR6825] - Nuclear Marker and ChIP Grade (ab129195) at 1/1000 dilution

Lane 1 : Wild-type HeLa cell lysate

Lane 2 : KDM1A knockout HeLa cell lysate

Lysates/proteins at 20 µg per lane.

Performed under reducing conditions.


Predicted band size: 92 kDa

Observed band size: 110 kDa

Lanes 1- 2: Merged signal (red and green). Green - ab129195 observed at 110 kDa. Red - Anti-GAPDH antibody [6C5] - Loading Control (**ab8245**) observed at 37 kDa.

ab129195 was shown to react with KDM1/LSD1 in wild-type HeLa cells in western blot. Loss of signal was observed when knockout cell line **ab265790** (knockout cell lysate **ab256965**) was used. Wild-

type HeLa and KDM1A knockout HeLa cell lysates were subjected to SDS-PAGE. Membrane was blocked for 1 hour at room temperature in 0.1% TBST with 3% non-fat dried milk. ab129195 and Anti-GAPDH antibody [6C5] - Loading Control ([ab8245](#)) overnight at 4°C at a 1 in 1000 dilution and a 1 in 20000 dilution respectively. Blots were developed with Goat anti-Rabbit IgG H&L (IRDye®800CW) preadsorbed ([ab216773](#)) and Goat anti-Mouse IgG H&L (IRDye®680RD) preadsorbed ([ab216776](#)) secondary antibodies at 1 in 20000 dilution for 1 hour at room temperature before imaging.




Immunoprecipitation - Anti-KDM1/LSD1 antibody
[EPR6825] - Nuclear Marker and ChIP Grade
(ab129195)

ab129195 (purified) at 1/20 immunoprecipitating KDM1/LSD1 in 10 µg Jurkat cell lysate (Lanes 1 and 2, observed at 110 kDa). Lane 3 - Rabbit monoclonal IgG ([ab172730](#)).


For western blotting, VeriBlot for IP Detection Reagent (HRP) ([ab131366](#)), was used for detection at 1/10,000 dilution.

Blocking/Dilution buffer and concentration: 5% NFDm/TBST.



ChIP - Anti-KDM1/LSD1 antibody [EPR6825] -
Nuclear Marker and ChIP Grade (ab129195)


Chromatin was prepared from HCT 116 cells according to the Abcam X-ChIP protocol. Cells were fixed with 1% formaldehyde for 10 minutes. The ChIP was performed with 25µg of chromatin, 5µg of ab129195 (red), and 20µl of protein A/G sepharose beads slurry (10µl of sepharose A beads + 10µl of sepharose G beads). 5µg of rabbit normal IgG was added to the beads control (grey). The immunoprecipitated DNA was quantified by real time PCR (Sybr green approach).



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-KDM1A/LSD1 antibody [EPR6825] - Nuclear Marker and ChIP Grade (ab129195)


Immunohistochemical staining of paraffin embedded rat kidney with purified ab129195 at a working dilution of 1/50. The secondary antibody used is **ab97051**, a goat anti-rabbit IgG (H&L) at a dilution of 1/500. The sample is counter-stained with hematoxylin. Antigen retrieval was performed using Tris-EDTA buffer, pH 9.0.

PBS was used instead of the primary antibody as the negative control, and is shown in the inset.




Immunocytochemistry/ Immunofluorescence - Anti-KDM1A/LSD1 antibody [EPR6825] - Nuclear Marker and ChIP Grade (ab129195)

ab129195 staining KDM1A/LSD1 in wild-type HAP1 cells (top panel) and KDM1A knockout HAP1 cells (bottom panel). The cells were fixed with 4% formaldehyde (10min), permeabilized with 0.1% Triton X-100 for 5 minutes and then blocked with 1% BSA/10% normal goat serum/0.3M glycine in 0.1% PBS-Tween for 1h. The cells were then incubated with ab129195 at 1µg/ml concentration and **ab195889** at 1/250 dilution (shown in pseudo colour red) overnight at +4°C, followed by a further incubation at room temperature for 1h with a goat secondary antibody to Rabbit IgG (Alexa Fluor® 488) (**ab150081**) at 2 µg/ml (shown in green). Nuclear DNA was labelled in blue with DAPI.



Immunocytochemistry/ Immunofluorescence - Anti-KDM1/LSD1 antibody [EPR6825] - Nuclear Marker and ChIP Grade (ab129195)

Immunofluorescence staining of HeLa cells with purified ab129195 at a working dilution of 1/100, counter-stained with DAPI. The secondary antibody was Alexa Fluor® 488 goat anti-rabbit ([ab150077](#)), used at a dilution of 1/1000. [ab7291](#), a mouse anti-tubulin antibody (1/1000), was used to stain tubulin along with [ab150120](#) (Alexa Fluor® 594 goat anti-mouse, 1/1000), shown in the top right hand panel. The cells were fixed in 4% PFA and permeabilized using 0.1% Triton X 100. The negative controls are shown in bottom middle and right hand panels - for negative control 1, purified ab129195 was used at a dilution of 1/500 followed by an Alexa Fluor® 594 goat anti-mouse antibody ([ab150120](#)) at a dilution of 1/500. For negative control 2, [ab7291](#) (mouse anti-tubulin) was used at a dilution of 1/500 followed by an Alexa Fluor® 488 goat anti-rabbit antibody ([ab150077](#)) at a dilution of 1/400.



Western blot - Anti-KDM1/LSD1 antibody [EPR6825] - Nuclear Marker and ChIP Grade (ab129195)

All lanes : ab129195 at 1/1000 dilution

Lane 1 : Mouse brain lysates

Lane 2 : Mouse heart lysates

Lane 3 : Mouse liver lysates

Lane 4 : Mouse spleen lysates

Lysates/proteins at 20 µg per lane.

Secondary


All lanes : Goat Anti-Rabbit IgG, (H+L), Peroxidase conjugated ([ab97051](#)) at 1/20000 dilution

Predicted band size: 92 kDa

Observed band size: 110 kDa

Exposure time: 180 seconds

Blocking/diluting buffer and concentration: 5% NFDM/TBST



All lanes : ab129195 at 1/1000 dilution

Lane 1 : Jurkat (Human T cell leukemia T lymphocyte) whole cell lysates

Lane 2 : HCT 116 (Human colorectal carcinoma epithelial cell) whole cell lysates

Lane 3 : NIH/3T3 (Mouse embryonic fibroblast) whole cell lysates

Lane 4 : C2C12 (Mouse myoblasts myoblast) whole cell lysates

Lane 5 : Mouse skeletal muscle lysates

Lane 6 : C6 (Rat glial tumor glial cell) whole cell lysates

Lysates/proteins at 20 µg per lane.

Secondary


All lanes : Goat Anti-Rabbit IgG, (H+L), Peroxidase conjugated at 1/20000 dilution

Predicted band size: 92 kDa

Exposure time: 20 seconds

Blocking/diluting buffer and concentration: 5% NFDM/TBST

Observed band: 110kd



Western blot - Anti-KDM1/LSD1 antibody [EPR6825] - Nuclear Marker and ChIP Grade (ab129195)

All lanes : Anti-KDM1/LSD1 antibody [EPR6825] - Nuclear Marker and ChIP Grade (ab129195) at 1/10000 dilution

Lane 1 : Wild-type HAP1 cell lysate

Lane 2 : KMD1 / LSD1 knockout HAP1 cell lysate

Lane 3 : HeLa cell lysate


Lane 4 : Jurkat cell lysate

Lysates/proteins at 20 µg per lane.

Predicted band size: 92 kDa

Lanes 1 - 4: Merged signal (red and green). Green - ab129195 observed at 110 kDa. Red - loading control, **ab8245**, observed at 37 kDa.

ab129195 was shown to specifically react with KMD1 / LSD1 in wild-type HAP1 cells. No band was observed when KMD1 / LSD1 knockout samples were used. Wild-type and KMD1 / LSD1 knockout samples were subjected to SDS-PAGE. ab129195 and **ab8245** (loading control to GAPDH) were both diluted 1/10,000 and incubated overnight at 4°C. Blots were developed with Goat anti-Rabbit IgG H&L (IRDye® 800CW) preadsorbed (**ab216773**) and Goat anti-Mouse IgG H&L (IRDye® 680RD) preadsorbed (**ab216776**) secondary antibodies at 1/10,000 dilution for 1 hour at room temperature before imaging.



Western blot - Anti-KDM1/LSD1 antibody [EPR6825] - Nuclear Marker and ChIP Grade (ab129195)

All lanes : Anti-KDM1/LSD1 antibody [EPR6825] - Nuclear Marker and ChIP Grade (ab129195) at 1/10000 dilution (purified)

Lane 1 : C6 whole cell lysate

Lane 2 : Raw 264.7 whole cell lysate

Lane 3 : PC-12 whole cell lysate

Lane 4 : NIH/3T3 whole cell lysate

Lysates/proteins at 10 µg per lane.


Secondary

All lanes : Goat Anti-Rabbit IgG H&L (HRP) (**ab97051**) at 1/20000 dilution

Predicted band size: 92 kDa


Observed band size: 110 kDa

Blocking/Dilution buffer: 5% NFDM/TBST.



Flow Cytometry (Intracellular) - Anti-KDM1/LSD1 antibody [EPR6825] - Nuclear Marker and ChIP Grade (ab129195)


Intracellular Flow Cytometry analysis of HeLa (human cervix adenocarcinoma) cells labeling KDM1/LSD1 with purified ab129195 at 1/20 dilution (10ug/mL) (red). Cells were fixed with 4% paraformaldehyde and permeabilised with 90% methanol. A Goat anti rabbit IgG (Alexa Fluor® 488) (1/2000 dilution) was used as the secondary antibody. Rabbit monoclonal IgG (Black) was used as the isotype control, cells without incubation with primary antibody and secondary antibody (Blue) were used as the unlabeled control.



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-KDM1/LSD1 antibody [EPR6825] - Nuclear Marker and ChIP Grade (ab129195)

Immunohistochemical staining of paraffin embedded mouse colon with purified ab129195 at a working dilution of 1/50. The secondary antibody used is **ab97051**, a goat anti-rabbit IgG (H&L) at a dilution of 1/500. The sample is counter-stained with hematoxylin. Antigen retrieval was performed using Tris-EDTA buffer, pH 9.0.

PBS was used instead of the primary antibody as the negative control, and is shown in the inset.



All lanes : Anti-KDM1/LSD1 antibody [EPR6825] - Nuclear Marker and ChIP Grade (ab129195) at 1/10000 dilution (purified)

Lane 1 : HEK293 whole cell lysate

Lane 2 : HeLa whole cell lysate

Lane 3 : Jurkat whole cell lysate

Lysates/proteins at 10 µg per lane.


Secondary

All lanes : Goat Anti-Rabbit IgG H&L (HRP) (**ab97051**) at 1/20000 dilution

Predicted band size: 92 kDa


Observed band size: 110 kDa

Blocking/Dilution buffer: 5% NFDM/TBST.



Immunohistochemical staining of paraffin embedded human stomach carcinoma with purified ab129195 at a working dilution of 1/50. The secondary antibody used is **ab97051**, a goat anti-rabbit IgG (H&L) at a dilution of 1/500. The sample is counter-stained with hematoxylin. Antigen retrieval was performed using Tris-EDTA buffer, pH 9.0.


PBS was used instead of the primary antibody as the negative control, and is shown in the inset.



Immunocytochemistry/ Immunofluorescence - Anti-KDM1/LSD1 antibody [EPR6825] - Nuclear Marker and ChIP Grade (ab129195)

This image is courtesy of an anonymous Abreview

Unpurified ab129195 staining KDM1/LSD1 in paraffin-embedded A549 lung cancer cells by ICC/IF (Immunocytochemistry/immunofluorescence). Cells were fixed using the HOPE technique and permeabilized with 0.05% Tween. Samples were incubated with primary antibody (1/100) for 45 minutes at 25°C. An Alexa Fluor®488-conjugated Donkey anti-mouse IgG polyclonal (1/200) was used as the secondary antibody.



Western blot - Anti-KDM1/LSD1 antibody [EPR6825] - Nuclear Marker and ChIP Grade (ab129195)

All lanes : Anti-KDM1/LSD1 antibody [EPR6825] - Nuclear Marker and ChIP Grade (ab129195) at 1/10000 dilution (unpurified)

Lane 1 : 293T cell lysate

Lane 2 : HeLa cell lysate

Lane 3 : Jurkat cell lysate

Lane 4 : PC3 cell lysate

Lane 5 : C6 cell lysate

Lane 6 : RAW264.7 cell lysate

Lane 7 : PC12 cell lysate


Lane 8 : NIH 3T3 cell lysate

Lysates/proteins at 10 µg per lane.

Secondary

All lanes : Goat anti-Rabbit HRP at 1/2000 dilution

Predicted band size: 92 kDa



Unpurified ab129195, at 1/100, staining KDM1/LSD1 in paraffin embedded human testis tissue by Immunohistochemistry.

Perform heat mediated antigen retrieval before commencing with IHC staining protocol.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-KDM1/LSD1 antibody [EPR6825] - Nuclear Marker and ChIP Grade (ab129195)

Why choose a recombinant antibody?



Research with confidence
Consistent and reproducible results



Long-term and scalable supply
Recombinant technology



Success from the first experiment
Confirmed specificity



Ethical standards compliant
Animal-free production

Anti-KDM1/LSD1 antibody [EPR6825] - Nuclear Marker and ChIP Grade (ab129195)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors