

Anti-KLRA1 antibody ab83946

3 图像

概述	
产品名称	Anti-KLRA1抗体
描述	兔多克隆抗体to KLRA1
宿主	Rabbit
经测试应用	适用于: ICC/IF, WB, IHC-P
种属反应性	与反应: Human
免疫原	13 amino acid synthetic peptide designed within the C terminal region of Human KLRA1
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

性能	
形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	pH: 7.2 Preservative: 0.02% Sodium azide Constituent: PBS
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用	
The Abpromise guarantee	Abpromise™ 承诺保证使用ab83946于以下的经测试应用
“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。	


应用	Ab评论	说明
ICC/IF		Use a concentration of 20 µg/ml.
WB		Use a concentration of 1 - 2 µg/ml. Predicted molecular weight: 25 kDa.
IHC-P		Use a concentration of 2.5 µg/ml.

靶标

相关性


KLRA1 is one of a cluster of lectin genes. It is likely to represent an evolutionarily remnant.

图片



ab83946 at 2.5µg/ml staining KLRA1 in human spleen by Immunohistochemistry.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-KLRA1 antibody (ab83946)



Lane 1 : Anti-KLRA1 antibody (ab83946) at 1 µg/ml

Lane 2 : Anti-KLRA1 antibody (ab83946) at 2 µg/ml

All lanes : 293 cell lysate

Lysates/proteins at 15 µg per lane.

Predicted band size: 25 kDa

Observed band size: 27 kDa

Western blot - Anti-KLRA1 antibody (ab83946)

Immunocytochemistry/ Immunofluorescence - Anti-KLRA1 antibody (ab83946)

Immunofluorescence of KLRA1 in Human Spleen cells using ab83946 at 20 ug/ml.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors