

# Anti-L-Alanine antibody ab9436

## 1 References

### 概述

产品名称	Anti-L-Alanine抗体
描述	兔多克隆抗体to L-Alanine
宿主	Rabbit
特异性	The antibody is calibrated against a spectrum of antigens to assure hapten selectivity. No measurable cross-reactivity (<1:1000) was detected against alanine in peptides or proteins. Fixed tissue cross-reactivity was tested with known targets at the recommended dilution. No measurable glutaraldehyde-fixed tissue cross-reactivity (<1:1000) was detected against D/L-aspartate, D/L-arginine, agmatine, guanidine, gamma-aminobutyrate, L-citrulline, L-cysteine, D/L-glutamate, D/L-glutamine, glutathione, glycine, L-lysine, L-ornithine, L-serine, taurine, L-threonine, L-tryptophan, L-tyrosine. Significant cross-reactivity (1:10) was detected against beta-alanine.
经测试应用	适用于: IHC-FoFr, ICC, IHC-Fr, Immunomicroscopy
种属反应性	与反应: Species independent
免疫原	Chemical/ Small Molecule corresponding to L-Alanine (Glutaraldehyde).
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at +4°C. Do Not Freeze.
纯度	IgG fraction
克隆	多克隆
同种型	IgG

### 应用

## The Abpromise guarantee

**Abpromise™** 承诺保证使用ab9436于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
IHC-FoFr		
ICC		
IHC-Fr		
Immunomicroscopy		

### 应用说明

This antibody is recommended for thin section post-embedding endogenous content mapping by Light and Electron microscopy immunostaining via High Performance Immunocytochemistry (HPI) (see "Protocols" tab above)

Recommended user dilution: 1/100  
True Dilution at User Dilution: 1/8000

High-performance thin section immunostaining using silver-intensified immunogold or fluorescence detection.

Enhanced detection is possible using streptavidin detection.

EM applications with 10-40 nm gold GAR IgGs or gold streptavidin.

OPTIMAL FIXATION: 0.1-2.5% glutaraldehyde, 1% formaldehyde using HPI protocol (The antisera targets the glutaraldehyde conjugate of the hapten).

MINIMAL GLUTARALDEHYDE: 0.1% using HPI

The hapten is osmium tolerant (deosmication required), therefore the antisera can be used on conventional post-embedding electron microscope immunostaining.

All procedures may be carried out at room temperature. Exact dilutions for all applications cannot be predicted, but it is unlikely that deviations from the calibrated levels will be needed. Dilutions are optimized for antigen detection over a 2 log unit range.

The product is optimized for HPI/EHPI with gold or fluorescence detection using etched plastic sections. Filter diluted reagents with 0.2 mm syringe filters before use on EM grids. Enzyme-linked visualizations can be used but will compress the signal dynamic range and are less sensitive.

Use with frozen or vibratome sections is possible but will not yield optimal images as IgGs penetrate aldehyde cross-linked tissue poorly and most amino acids are present at such high levels that prozone effects occur. Use in whole mounts is not recommended for similar reasons.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet

- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors