

### Anti-L3MBTL3 antibody ab68117

★★★★★ [2 Abreviews](#) [2 图像](#)

#### 概述

产品名称	Anti-L3MBTL3抗体
描述	小鼠多克隆抗体to L3MBTL3
宿主	Mouse
经测试应用	适用于: WB, ICC/IF
种属反应性	与反应: Human
免疫原	Recombinant full length protein within Human L3MBTL3. The exact immunogen sequence used to generate this antibody is proprietary information. If additional detail on the immunogen is needed to determine the suitability of the antibody for your needs, please <b><a href="#">contact</a></b> our Scientific Support team to discuss your requirements.

#### 常规说明

The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.

If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As

#### 性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	pH: 7.40 Constituent: 100% PBS
纯度	Protein G purified
克隆	多克隆
同种型	IgG

#### 应用

The Abpromise guarantee


Abpromise™承诺保证使用ab68117于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB	★★★★★ (1)	1/500 - 1/1000. Detects a band of approximately 100 kDa (predicted molecular weight: 88 kDa).
ICC/IF		Use a concentration of 10 µg/ml.

靶标	
功能	Putative Polycomb group (PcG) protein. PcG proteins maintain the transcriptionally repressive state of genes, probably via a modification of chromatin, rendering it heritably changed in its expressibility. Required for normal maturation of myeloid progenitor cells.
序列相似性	Contains 3 MBT repeats. Contains 1 SAM (sterile alpha motif) domain.
细胞定位	Nucleus.

图片



Western blot - Anti-L3MBTL3 antibody (ab68117)

**All lanes** : Anti-L3MBTL3 antibody (ab68117) at 1/500 dilution

**Lane 1** : L3MBTL3 transfected 293T cell lysate

**Lane 2** : Non transfected 293T cell lysate


Lysates/proteins at 25 µg per lane.

**Secondary**

**All lanes** : Goat Anti-Mouse IgG (H&L)-HRP Conjugate secondary antibody at 1/2500 dilution

**Predicted band size:** 88 kDa

**Observed band size:** 100 kDa



Immunofluorescence analysis of HeLa cells using purified ab68117 at 10 μg/ml followed by fluorescein conjugated secondary antibody. Cells were fixed with 4% paraformaldehyde and permeabilised with 0.1% Triton X-100.

Immunocytochemistry/ Immunofluorescence - Anti-L3MBTL3 antibody (ab68117)

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors