

# Anti-Leptospira biflexa antibody ab20110

## 概述

产品名称	Anti-Leptospira biflexa抗体
描述	兔多克隆抗体to Leptospira biflexa
宿主	Rabbit
特异性	This antibody is specific for all antigens of Leptospira biflexa.
经测试应用	适用于: ICC/IF
种属反应性	与反应: Leptospira biflexa
免疫原	Tissue, cells or virus corresponding to Leptospira biflexa. Leptospira biflexa, whole organism
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

## 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: 0.1% Sodium azide Constituent: 0.0268% PBS
纯度	Protein A purified
纯化说明	>95% pure.
克隆	多克隆
同种型	IgG

## 应用

The Abpromise guarantee

**Abpromise™** 承诺保证使用ab20110于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
ICC/IF		Use at an assay dependent concentration.

## 靶标

相关性	Leptospirosis, a contagious disease affecting both animals and humans and spread by infection with a bacterial pathogen called Leptospira, may result in chronic liver and kidney disease and fatality in the dog. Leptospira is a flexible, spiral-shaped, Gram-negative spirochete with internal flagella. It thrives in water and appears long and helical with a characteristic hook on one or both ends. These organisms are divided into two species, Leptospira biflexa and Leptospira interrogans.
-----	--

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors