

# Anti-LRP5L antibody ab122013

## 2 图像

### 概述

产品名称	Anti-LRP5L 抗体
描述	兔多克隆抗体 to LRP5L
宿主	Rabbit
经测试应用	适用于: ICC/IF, IHC-P
种属反应性	与反应: Human
免疫原	amino acids 198-251 (VIIDQLPDLM GLKAVNVDKV VGTNPHADRN GGAATCASSR PTQPGLAAPS RAWN) of Human LRP5L (A4QPB2).  <a href="#">Run BLAST with ExPASy</a>  <a href="#">Run BLAST with NCBI</a>
阳性对照	ICC: U-251 MG cells; IHC-P: Human skeletal muscle tissue.
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

### 性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
存储溶液	pH: 7.20 Preservative: 0.02% Sodium azide Constituents: 40% Glycerol (glycerin, glycerine), 59% PBS
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

### 应用

## The Abpromise guarantee

**Abpromise™** 承诺保证使用ab122013于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。


应用	Ab评论	说明
ICC/IF		Use a concentration of 0.25 - 2 µg/ml.
IHC-P		1/20 - 1/50. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.

## 靶标

### 序列相似性

Contains 5 LDL-receptor class B repeats.

## 图片



Immunocytochemistry/ Immunofluorescence - Anti-LRP5L antibody (ab122013)

Immunofluorescent staining of Human cell line U-251 MG shows positivity in nucleus but not nucleoli. Recommended concentration of ab122013 1-4 µg/ml. Cells treated with PFA/Triton X-100.



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-LRP5L antibody (ab122013)

ab122013, at a 1/10 dilution, staining Human LRP5L in skeletal muscle, using Immunohistochemistry, Formalin/PFA-fixed paraffin-embedded tissue. Moderate cytoplasmic positivity in myocytes shown.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors