

# Anti-Lysozyme antibody ab34799

★★★★☆ [1 Abreviews](#) [2 References](#) [1 图像](#)

### 概述

---

产品名称	Anti-Lysozyme抗体
描述	兔多克隆抗体to Lysozyme
宿主	Rabbit
特异性	ab34799 recognises purified and partially purified Lysozyme.
经测试应用	适用于: WB
种属反应性	与反应: Chicken
免疫原	Full length native protein (purified) corresponding to Chicken Lysozyme. Purified Lysozyme from chicken egg white. Database link: <a href="#">P00698</a>
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

### 性能

---

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term.
存储溶液	Preservative: 0.01% Sodium azide Constituents: 0.88% Sodium chloride, 0.42% Potassium phosphate
纯化说明	Purified from monospecific antiserum by a delipidation and defibrination.
克隆	多克隆
同种型	IgG

### 应用

---

## The Abpromise guarantee

**Abpromise™** 承诺保证使用 ab34799 于以下的经测试应用


“应用说明”部分下显示的仅为推荐的起始稀释度；实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
WB	★★★★☆ (1)	1/500 - 1/5000. Predicted molecular weight: 17 kDa.

## 靶标

功能	Lysozymes have primarily a bacteriolytic function; those in tissues and body fluids are associated with the monocyte-macrophage system and enhance the activity of immunoagents.
疾病相关	Amyloidosis 8
序列相似性	Belongs to the glycosyl hydrolase 22 family.
细胞定位	Secreted.

## 图片



Anti-Lysozyme antibody (ab34799) at 1/1000 dilution + Lysozyme at 0.01 µg

### Secondary

Goat anti-rabbit DyLight 649 at 1/20000 dilution

**Predicted band size:** 17 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors