

### Anti-MAD1L1/MAD1 antibody [9B10] ab5783

★★★★★ [2 Abreviews](#) [5 References](#) [1 图像](#)

#### 概述

产品名称	Anti-MAD1L1/MAD1抗体[9B10]
描述	小鼠单克隆抗体[9B10] to MAD1L1/MAD1
宿主	Mouse
经测试应用	<b>适用于:</b> Flow Cyt <b>不适用于:</b> WB
种属反应性	<b>与反应:</b> Human
免疫原	Recombinant full length protein corresponding to Human MAD1L1/MAD1. Database link: <a href="#">Q9Y6D9</a> , <a href="#">Q13586</a>
常规说明	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

#### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: 0.1% Sodium azide Constituent: PBS
纯度	Protein A purified
克隆	单克隆
克隆编号	9B10
同种型	IgG2b

#### 应用

The Abpromise guarantee

Abpromise™承诺保证使用ab5783于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。

应用	Ab评论	说明
Flow Cyt		Use 0.1µg for 10 <sup>6</sup> cells. <b>ab170192</b> - Mouse monoclonal IgG2b, is suitable for use as an isotype control with this antibody.


应用说明

Is unsuitable for WB.

靶标

功能	Component of the spindle-assembly checkpoint that prevents the onset of anaphase until all chromosomes are properly aligned at the metaphase plate. May recruit MAD2L1 to unattached kinetochores. Has a role in the correct positioning of the septum. Required for anchoring MAD2L1 to the nuclear periphery. Binds to the TERT promoter and represses telomerase expression, possibly by interfering with MYC binding.
组织特异性	Expressed weakly at G0/G1 and highly at late S and G2/M phase.
疾病相关	Defects in MAD1L1 are involved in the development and/or progression of various types of cancer.
序列相似性	Belongs to the MAD1 family.
翻译后修饰	Phosphorylated; by BUB1. Become hyperphosphorylated in late S through M phases or after mitotic spindle damage.
细胞定位	Nucleus. Chromosome > centromere > kinetochore. Cytoplasm > cytoskeleton > microtubule organizing center > centrosome. Cytoplasm > cytoskeleton > spindle. From the beginning to the end of mitosis, it is seen to move from a diffusely nuclear distribution to the centrosome, to the spindle midzone and finally to the midbody. Colocalizes with NEK2 at the kinetochore.

图片



Overlay histogram showing HeLa cells stained with ab5783 (red line). The cells were fixed with 80% methanol (5 min) and then permeabilized with 0.1% PBS-Tween for 20 min. The cells were then incubated in 1x PBS / 10% normal goat serum / 0.3M glycine to block non-specific protein-protein interactions followed by the antibody (ab5783, 0.1µg/1x10<sup>6</sup> cells) for 30 min at 22°C. The secondary antibody used was Alexa Fluor® 488 goat anti-mouse IgG (H+L) (**ab150113**) at 1/2000 dilution for 30 min at 22°C. Isotype control antibody (black line) was mouse IgG2b [PLPV219] (**ab91366**, 1µg/1x10<sup>6</sup> cells) used under the same conditions. Unlabelled sample (blue line) was also used as a control. Acquisition of >5,000 events were collected using a 20mW Argon ion laser (488nm) and 525/30 bandpass filter.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors